

## Product datasheet for **MR226176L3V**

### **IL12b (NM\_008352) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	IL12b (NM_008352) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	IL12b
Synonyms:	IL-12b; IL-12p40; IL12p40; p40
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008352
ORF Size:	1005 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226176).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_008352.2</a> , <a href="#">NP_032378.1</a>
RefSeq Size:	1951 bp
RefSeq ORF:	1008 bp
Locus ID:	16160
Cytogenetics:	11 25.94 cM
Gene Summary:	This gene encodes the beta subunit p40 of the Interleukin 12 (IL-12) family of cytokines. Members of the IL-12 family form heterodimers consisting of heavy and light subunits linked by disulfide bonds. The product of this gene, p40, is a subunit of interleukins IL-12 and IL-23. [provided by RefSeq, Dec 2014]


[View online »](#)