

Product datasheet for **MR226183**

Zbtb16 (NM_001033324) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zbtb16 (NM_001033324) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Zbtb16
Synonyms:	AI467657; lu; PLZF; Zfp145
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR226183 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGGATCTGACAAAGATGGGGATGATCCAGCTGCAGAACCTAGCCACCCACGGGGCTGCTGTGCAAGG
CCAATCAGATGCGACTGGCTGGGACTTTGTGCGATGTGGTCATCATGGTGACAGCCAGGAGTTCCATGC
CCACCGGACGGTCTAGCCTGCACCAGCAAGATGTTTGAGATCCTCTTCCACCGAAAACAGCCAGCACTAT
ACTCTGGACTTCTCTCGCCAAAAACCTTCCAGCAGATCCTGGAGTACGCCTACACGGCCACACTGCAAG
CCAAGGCGGAGGACCTGGATGACCTGCTGTATGCAGCTGAGATTTAGAGATCGAATACCTGGAGGAGCA
GTGCTGAAGATCCTGGAGACCACCCAGGCATCTGATGACAATGACACAGAGGCCACCATGGCTGACGGT
GGGGGCGAAGAAGAAGAGGACCGTAAGGCTCGATACCTCAAGAACATCTTTATCTCGAAGCATTCCAGCG
AGGAGAGTGGCTACGCCAGTGTGGCTGGACAGAGCCTCCCTGGGCCATGGTGGACCAGAGCCCCCTCAGT
CTCCACCTCTTTCGGTCTCTCAGCCATGAGTCCTACCAAGGCAGCTGTGGACAGCTTGATGAGTATAGGA
CAGTCACTCCTGCAAGGAACTTTCAGCCACCTGCAGGGCCTGAGGAGCCACACTGGCAGGGGGTGGGC
GACACCCTGGGGTGGCTGAGGTGAAGATGGAAATGATGCAGGTAGATGAAGCCCTGGCCAGGACAGCCC
TGGGGCAGCTGAGTCTAGCATCTCAGGAGGGATGGGGGACAAGCTTGAAGAGAGGAGCAAAGAGGGGCT
GGGACCCCGACTCGAGGCAGCGTCATCACTAGTGCCAGAGAGCTGCATTATGGGAGAGAGGAGTGGTG
AGCAGCTCTCGCTCCTGTTGAAGCTGGCCAGGGACCCCTGGGGCGCAGGAGCCCTGGCACCCCCAGT
GGAGAAGCATTGGGTATCTACTCGGTGCTGCCAACCACAAGGCCGATGCTGTGTTGAGCATGCCGTCT
TCAGTGACGTCTGGTCTCCATGTGCAACCTGCCCTGGCAGTCTCCATGGACTTCAGCACCTACGGGGGTC
TGCTGCCCTCAGGGCTTCATCCAGAGGGAGCTGTTGAGCAAGCTGGGGGAGCTGGCTGTGGCATGAAGGC
TGAGAGCCGCCCTCTGGGGGAGCAGTGCAGTGTGTGTGGGGTTCGAGCTTCCGGACAACGAGGCAGTGGAG
CAGCACAGGAACTGCACAGTGGGATGAAAACATACGGGTGTGAACTCTGCGGAAAACGGTTCCTGGACA
GTTTGGGACTGAGAATGCATTTACTGGCTCATTACGCGGGTGCCAAAGCCTTTGTGTGTGATCAGTGTGG
TGCCAGTTCTCAAAGGAGGATGCCCTGGAGACACACAGACAGACCCATACTGGCACGGACATGGCTGTC
TTCTGTCTGCTGTGGGAAACGCTTTCAGGCACAAAGCGCACTCCAGCAGCACATGGAGGTCCACGCAG
GCGTGCGCAGCTATATTTGCAGTGAGTGAACCGCACCTCCCCAGCCACACGGCTCTCAAGCGCCACCT
TCGCTCACATACAGGTGACCACCATATGAGTGTGAGTCTGCGGCAGCTGCTTCCGGGATGAGAGCACA
CTCAAGAGCCACAAGCGCATCCACACAGGGGAGAAGCCCTATGAGTGTAAACGGCTGTGGCAAGAAGTTCA
GCCTCAAGCACCAAGTGGAGACGCACTACAGGGTTCACACAGGTGAGAAGCCCTTTGAGTGCAAACCTCTG
CCACCAGCGCTCCCGAGACTACTCGGCCATGATCAAGCACCTGAGAACACACAACGGTGCCTCGCCATAC
CAGTGTACCATCTGCACGGAATACTGCCCCAGCCTCTCCTCCATGCAGAAAACACATGAAGGGCCACAAGC
CGGAGGAGATCCCACCTGACTGGAGGATAGAGAAGACGTACCTCTACCTGTGTTATGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR226183 protein sequence
Red=Cloning site Green=Tags(s)

```
MDLTKMGMQLQNPSHTGLLCKANQMRLAGTLCDDVIMVDSQEFHAHRTVLACTSKMFEILFHRNSQHY
TLDLFLSPKTFQQILEYAYATLQAKAEDLDDL YAAEILEIEYLEEQCLKILETIQASDDNDTEATMADG
GGEEEEERKARYLKNIFI SKHSSEESGYASVAGQSLPGPMVDQSPSVSTSFGLSAMSPTKAAVDSLMSIG
QSLQLQGTLPAPGPEEPTLAGGGRHPGVAEVKMEMMQVDEAPGQDSPGAAESSISGGMGDKLEERSKEGP
GTPTRGSVITSARELHYGREESGEQLSPPVEAGQPPGRQEPLAPPVEKHLGIYSVLPNHKADAVLSMPS
SVTSGLVHPALAVSMDFSTYGGLLPQGFIQRELF SKLGELAVGMKAESRPLGEQCSVCVGPVLPDNEAVE
QHRKLSHGKMTYGCCELCKRFLDSLRLRMHLLAHSAGAKAFVCDQCGAQF SKEDALETHTQHTGTDMAV
FCLLCGKRFQAQSALQQHMEVHAGVRSYICSECNRTFPSHTALKRHLRSHTGDHPYCEFCGSCFRDEST
LKSHKRIHTGEKPYECNGCGKFKSLKHQLETHYRVHTGEKPFECKLCHQSRDYSAMIKHLRTHNGASPY
QCTICTEYCPSSLSSMQHKMGKHPPEEIPPDWRIEKTLYLYCYV
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001033324

ORF Size: 2022 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001033324.1](#), [NM_001033324.2](#), [NP_001028496.1](#)

RefSeq Size: 5114 bp


RefSeq ORF: 2022 bp

Locus ID: 235320

Cytogenetics: 9 26.47 cM

MW: 74.2 kDa

Product images:



Circular map for MR226183