

Product datasheet for **MR226185L1V**

Ccne1 (NM_007633) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ccne1 (NM_007633) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ccne1
Synonyms:	AW538188; CycE1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_007633
ORF Size:	1224 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226185).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007633.2 , NP_031659.2
RefSeq Size:	2022 bp
RefSeq ORF:	1227 bp
Locus ID:	12447
UniProt ID:	Q61457
Cytogenetics:	7 25.35 cM
Gene Summary:	Essential for the control of the cell cycle at the G1/S (start) transition.[UniProtKB/Swiss-Prot Function]



[View online »](#)