



Product datasheet for MR226204L4

Epha6 (NM_007938) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Epha6 (NM_007938) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Epha6
Synonyms:	Ehk2; Hek12; m-ehk2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226204).
Restriction Sites:	Sgfl-NotI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_007938
ORF Size:	3390 bp



[View online »](#)

