

Product datasheet for **MR226205L4V**

Hes5 (NM_010419) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hes5 (NM_010419) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hes5
Synonyms:	bHLHb38
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010419
ORF Size:	504 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226205).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010419.4 , NP_034549.1
RefSeq Size:	1273 bp
RefSeq ORF:	504 bp
Locus ID:	15208
UniProt ID:	P70120
Cytogenetics:	4 86.01 cM
Gene Summary:	Transcriptional repressor of genes that require a bHLH protein for their transcription. Plays an important role as neurogenesis negative regulator.[UniProtKB/Swiss-Prot Function]



[View online »](#)