

Product datasheet for **MR226242**

Nrl (NM_001136074) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nrl (NM_001136074) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nrl
Synonyms:	D14H14S46E
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR226242 representing NM_001136074 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGCTTTCCCTCCAGTCCCTTGGCTATGGAATATGTTAATGACTTTGATTTGATGAAGTTCGAAATAA
AGCGGGAGCCTTCTGAGGGCCGATCTGGAGTCCCCACAGCCTCGTTGGGCTCCACACCATACAGCTCGGT
GCCTCCTTCACCACCTTCAGTGAGCCAGGCATGGTGGTGGTGGCGAGGCCCTAGGCCAGGCCTGGAG
GAGCTATATTGGCTGGCCACCCTGCAGCAGCAGCTGGGTGGATGAGGTTCTCGGGCTGAGTCCCGACG
AAGCTGTGGAGCTGCTGCAGAACAGGGTCCTGTCTCTATGGAAGGGCCTTTGGCTACTATTAGGGAG
CCCGGGAGAGACAGGAGCCAGCATGTCCAGCTGCCCAGAGATTTTCGGACGCCGCGCTGGTCTCGATG
TCTGTGCGCGAGTTGAACCGGCAGCTGCGGGGCTGCGGGCGCGACGAGGCTCTGCGGCTGAAGCAGAGGC
GCCGCACGCTCAAGAACCGCGGTACGCTCAGGCTTGTGCTCCAAGCGGCTGCAACAGCGGCGCGGGCT
GGAGGCAGAGCGCGCTCGCTGGCCGCCAGCTGGATGCCCTGCGGGCCGAGGTGGCAGCGCTGGCTCGC
GAGCGCACCTCTATAAGGCCGCTGTGACCGGCTGACCTCCGGCGGCCCTGGGTCCGACGACCACACAC
ACCTCTTCCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR226242 representing NM_001136074
 Red=Cloning site Green=Tags(s)


MAFPSPPLAMEYVNDLDMKFEIKREPSEGRSGVPTASLGSTPYSSVPPSPPTFSEPGMVGGGEAPRPGLE
 ELYWLATLQQQLGSDEVLGLSPDEAVELLQNQGPVSMEGPLGYSGSPGETGAQHVLPERFSDAALVSM
 SVRELNRQLRGCGRDEALRLKQRRRTLKNRGYAQACRSKRLQRRGLEAERARLAAQLDALRAEVARLAR
 ERDLYKARCDRLTSGGPGSDDHHLFL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9075_g12.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001136074

ORF Size: 711 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001136074.2](#), [NP_001129546.1](#)

RefSeq Size: 2484 bp

RefSeq ORF: 714 bp

Locus ID: 18185


UniProt ID: [P54846](#)

Cytogenetics: 14 28.19 cM

MW: 26.1 kDa

Gene Summary: This gene encodes a member of the basic leucine zipper domain family of transcription factors. The encoded protein is preferentially expressed in the retina and is necessary for rod photoreceptor development. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2012]

Product images:



Circular map for MR226242