

Product datasheet for **MR226242L4V**

Nrl (NM_001136074) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nrl (NM_001136074) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nrl
Synonyms:	D14H14S46E
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001136074
ORF Size:	711 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226242).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001136074.2 , NP_001129546.1
RefSeq Size:	2484 bp
RefSeq ORF:	714 bp
Locus ID:	18185
UniProt ID:	P54846
Cytogenetics:	14 28.19 cM


[View online »](#)

Gene Summary:

This gene encodes a member of the basic leucine zipper domain family of transcription factors. The encoded protein is preferentially expressed in the retina and is necessary for rod photoreceptor development. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2012]