

Product datasheet for MR226243L3

Slurp1 (NM_020519) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slurp1 (NM_020519) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Slurp1
Synonyms:	1110021N19Rik; AI415082; ARS; ArsB; Slurp-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226243).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_020519
ORF Size:	333 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_020519.1
RefSeq Size:	525 bp
RefSeq ORF:	333 bp
Locus ID:	57277
UniProt ID:	Q9Z0K7
Cytogenetics:	15 D3
Gene Summary:	<p>Has an antitumor activity. Was found to be a marker of late differentiation of the skin. Implicated in maintaining the physiological and structural integrity of the keratinocyte layers of the skin. In vitro down-regulates keratinocyte proliferation; the function may involve the proposed role as modulator of nicotinic acetylcholine receptors (nAChRs) activity. In vitro inhibits alpha-7-dependent nAChR currents in an allosteric manner (By similarity). In T cells may be involved in regulation of intracellular Ca(2+) signaling (PubMed:17286989). Seems to have a immunomodulatory function in the cornea. The function may implicate a possible role as a scavenger receptor for PLAU thereby blocking PLAU-dependent functions of PLAU such as in cell migration and proliferation (PubMed:23139280, PubMed:25168896).</p> <p>[UniProtKB/Swiss-Prot Function]</p>

Product images:



Circular map for MR226243L3