

Product datasheet for **MR226276L1V**

Efna4 (NM_007910) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Efna4 (NM_007910) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Efna4
Synonyms:	EFL-4; Epl4; LERK-4
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_007910
ORF Size:	618 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226276).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007910.2 , NP_031936.2
RefSeq Size:	1596 bp
RefSeq ORF:	621 bp
Locus ID:	13639
UniProt ID:	O08542
Cytogenetics:	3 39.07 cM



[View online »](#)

Gene Summary:

Cell surface GPI-bound ligand for Eph receptors, a family of receptor tyrosine kinases which are crucial for migration, repulsion and adhesion during neuronal, vascular and epithelial development. Binds promiscuously Eph receptors residing on adjacent cells, leading to contact-dependent bidirectional signaling into neighboring cells. May play a role in the interaction between activated B-lymphocytes and dendritic cells in tonsils (By similarity). [UniProtKB/Swiss-Prot Function]