

Product datasheet for MR226280L3V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Nr4a2 (NM_013613) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: Nr4a2 (NM 013613) Mouse Tagged ORF Clone Lentiviral Particle

Symbol: Nr4a2

Synonyms: HZF-3; NOT; Nurr1; RNR-1; TINOR; TINUR

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Tag: Myc-DDK
ACCN: NM 013613

ORF Size: 1794 bp

ORF Nucleotide

Sequence:

The ORF insert of this clone is exactly the same as(MR226280).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through

naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: <u>NM 013613.2</u>, <u>NP 038641.1</u>

 RefSeq Size:
 3172 bp

 RefSeq ORF:
 1797 bp

 Locus ID:
 18227

 UniProt ID:
 Q06219

Cytogenetics: 2 31.66 cM







Gene Summary:

Transcriptional regulator which is important for the differentiation and maintenance of mesodiencephalic dopaminergic (mdDA) neurons during development. It is crucial for expression of a set of genes such as SLC6A3, SLC18A2, TH and DRD2 which are essential for development of mdDA neurons.[UniProtKB/Swiss-Prot Function]