

Product datasheet for **MR226300L4V**

Musk (NM_010944) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Musk (NM_010944) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Musk |
| Synonyms: | Mdk4; Mlk; Ns; Nsk1; Nsk2; Nsk3 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_010944 |
| ORF Size: | 2604 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR226300). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_010944.2 , NP_035074.2 |
| RefSeq Size: | 3354 bp |
| RefSeq ORF: | 2607 bp |
| Locus ID: | 18198 |
| UniProt ID: | Q61006 |
| Cytogenetics: | 4 31.87 cM |



[View online »](#)

Gene Summary:

This gene encodes a member of the protein tyrosine kinase family. The encoded protein is a type 1 receptor-like protein located in muscle membrane that is activated by the heparan sulfate proteoglycan agrin released by nerve cells. The encoded protein activates signaling cascades responsible for multiple aspects of motor neuron and muscle development, including organization of the postsynaptic membrane, synaptic gene transcription, patterning of skeletal muscle, anchoring of acetylcholinesterase, and guidance of motor axons. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]