

Product datasheet for **MR226322**

Klf9 (NM_010638) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Klf9 (NM_010638) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Klf9
Synonyms: 2310051E17Rik; AA589643; BTEB-1; Bteb1; ENSMUSG00000055440; Gm9971
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR226322 representing NM_010638
Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGGCGGCCGGAATTCGTGCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCCGCGGCCGCTACATGGACTTCGTGGCTGCCAGTGTCTGGTTTCCATCTCCAACCGCGCCCGC
TGCCGGAGCACGGGGCGCTCCGGAAGCCGAGCGGCTGCGACTACCTGAGCGGAGGTGACCAAGGAACA
CGGTGACCCGGGGACACCTGGAAGGATTATTGCACGCTGGTCACTATCGCCAAGAGCTTGTGGACCTC
AACAAATACCGACCCATCCAGACCCCTCGGTGTGCAGCGACAGTCTGGAGAGTCCCGATGAGGATATAG
GATCCGACAGCGAGTACCACCGAATCTGGTTCGAGTCTTCCACAGCCGGAGGAGAGACAGGATTC
TGGCAGCGCGCCAGCCACTCTCCCTCCTCCACTCTGGAGTGGCTTCAAGGGGAAACACGCCTCCGAA
AAGAGGCACAAGTGCCCTACAGTGGCTGTGGAAAGTCTATGGAAAATCCTCCCATCTTAAGCCATT
ACAGAGTGCATACAGGTGAACGGCCCTTCCCTGCACGTGGCCAGACTGCCTTAAAAAGTTCTCGCGCTC
GGATGAGCTGACCCGCACTACCGGACCCACACTGGGAAAAGCAGTTCCGTTGCCCACTGTGTGAGAAG
AGATTCATGAGGAGTGACCATCTACCAAGCATGCCCGGCGTACACCGAGTTCATCCCAGCATGATCA
AGAGATCAAAAAGGCTTTGCCAGCCCTTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_010638

ORF Size: 732 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_010638.5](#)

RefSeq Size: 3263 bp

RefSeq ORF: 735 bp

Locus ID: 16601
 Cytogenetics: 19 B
 MW: 27.6 kDa

Product images:



Circular map for MR226322