

## Product datasheet for **MR226324**

### Sh3kbp1 (NM\_021389) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Sh3kbp1 (NM_021389) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Sh3kbp1
Synonyms:	1200007H22Rik; 1700125L08Rik; 5830464D22Rik; AI447724; Cin85; IN85; Ruk; Seta
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide Sequence:**

>MR226324 representing NM\_021389  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGGAGCTTTCTGCAGCCAAAGCCCAAGCCCCACAGACTTGCCTGAGAGTGAAATAAAGAAAGATATGA  
 AGAAAGATCTTCTCTAACAAAGCTCCAGAGAAACCCATGCACGATGTATCCAGTGGAAACGCTCTGCT  
 TTCTTCTGAAACCATTTTAAAGACCAATAAACGAGGAGAACGACGGAGACGCCGATGCCAGTGGCATT  
 AGCTACCTGCCCCAGAATGACGATGAGCTGGAGCTGAAAGTCGGGGACATCATAGAGTTGTAGGAGAGG  
 TAGAGGAAGGATGGTGGGAAGGTGTGCTCAATGGGAAGACTGGAATGTTTCCTTCCAACCTCATCAAGGA  
 GCTGTCAGGGGAGTCGGATGAGCTTGGCATCTCCAGGATGAGCAGCTTTCGAAGTCAAGCTTAAGAGAA  
 ACCACAGGCTCCGAGAGTGATGGAGGCGACTCCAGCAGTACCAAGTCCGAAGGTGCCAACGGGACCATGG  
 CAACTGCAGCCATCCAGCCCAAGAAAGTAAAGGAGTGGGGTTTGAGACATTTTCAAAGACAAGCCAAAT  
 TAAACTAAGACCAAGGTCGATTGAAGTGGAAAATGACTTCTGCCCGTCGAAAAGACAATGGGAAGAAG  
 TTACCTCCAGCTACATCAACTCCAGACCCATCGAAGACGGAGATGGACAGCAGGACAAAGACCAAGGATT  
 ACTGCAAAGTAATATTTCCATATGAGGCACAGAATGATGATGAATTGACAAATCAAGGAAGGAGATATAGT  
 GACTCTCATCAATAAGGACTGCATCGACGTAGGCTGGTGGGAAGGAGAATCAATGGCAGACGAGGCGTG  
 TTTCTGACAACCTTCGTGAAGTACTTCCATCGGACTTTGACAAGGAAGGGAATAGACCCAAGAAACCGC  
 CACCTCCATCAGTCTCTGTCGTCAAACAAGGGGACGTTACCACAGAGAGAAAAACATGAAATTAAGAAAT  
 ACCTCTGAAAGACCAGAAACCTTCCAAACAGAACAGAGAGAAAGAAAGACCAGAGAGAGAGCCAAAA  
 CTGGATTTACAGAAGCCTTCAGTTCCTGCCATCCCACAAAAAGCCTCGGCCACCTAAGACCAATTCCTC  
 TGAACAGACTGGTGCATTGCCCAAGAAGACTGAGCGACAGTGGGCCACTGACACACACCAAGGGG  
 TGACAGTCCAAAGATTGACTTGGCAGGACGCGCACTCTCTGGCATCTGGATAAGGACCTCTCAGACCGC  
 AGCAATGATATTGACTTAGAAGTTTTGACTCTGTGATATCATCTACTGAGAACTCAGCCACCCAAACCA  
 CAAGCAGACCAAAAGCTACAGGAAGACGGCCTCCATCCCAGTCACTCACATCTTCTTCCCTTTCAAGCCC  
 CGATATCTTTGACTCCCCAAGTCCCGAAGAAGATAAAGAGGAACACATTTGCTTGGCAGACAGGAATA  
 GACGTGTCAAAGAAGACTTCCAAGACTGTACCATATCCCAAGTGTGATAACAAGACATCCCTGCCAC  
 CCAAGCCAGGGACCATGGCAGCTGCCAGCAGTGGGCCAGCTTCTCTCTTTCAGTGGCATCTCACCCAT  
 GTCATCTCTTTGGAACAGCTGGACAGAGAGCCAGTTCTCCATCTCTGTTGAGCAGAGAGGAAAGCCA  
 AAGATGGAGCCAGCAGTGAAGGACCAGCAGAAACGTGAGATTAAGCAGTACTGTGCAAGTCCGTGAGCTGAGGA  
 CCATCATTGAGACCATGAAGGACCAGCAGAAACGTGAGATTAAGCAGTACTGTGCAAGTCCGTGAGCTGAGGA  
 GAAAAAGATCCGGCTCCGGTTGCAGATGGAAGTGAACGACATAAAGAAAGCTCTTCAATCAAAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**


>MR226324 representing NM\_021389  
 Red=Cloning site Green=Tags(s)

MELSAKAPSPTDLPESEIKKDMKDLLSNKAPEKPMHDVSSGNALLSSETILRTNKRGERRRRRRCQVAF  
 SYLPQNDDLELELVGDIIEVVGEVEEGWEGVLNGKTMFSPNFIKELSGESDELGISQDEQLSKSSLRE  
 TTGSESDGGDSSSTKSEGANMTATAAIQPKKVKGVGFDIFKDKPIKLRPRRSIEVENDFLPVEKTIKGG  
 LPPATSTPDPKTEMDSRTKTKDYCKVIFPYEAQNDELTIKEGDIVTLINKDCIDVGGWEGELNGRRGV  
 FPDNFKVLLPSDFDKEGNRPKPPPPSAPVVKQGAGTTERKHEIKKIPPERPETLPNRTEEKERPEREPK  
 LDLQKPSVPAIPPKPRPKTNSLNRPGALPPRRPERPVGPLTHTRGDSPKIDLAGSALSGILDKDLSDR  
 SNDIDLEGFDSVISSTEKLSHPTTSRPKATGRRPPSQSLTSSSLSPDIFDSSPSPEEDKEEHI SLAHRGI  
 DVSKKTSKTVTISQVSDNKTSLPPKPGTMAAASSGPASLSSVASSPMSSSLGTAGQRASSPSLFFSTEGKP  
 KMEPAVSSQAAIEELKMQVRELRTIETMKDQKREIKQLLSELDEEKIRLRLQMEVNDIKKALQSK


**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_021389

ORF Size: 1884 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_021389.6</a> , <a href="#">NP_067364.2</a>
<b>RefSeq Size:</b>	4730 bp
<b>RefSeq ORF:</b>	1887 bp
<b>Locus ID:</b>	58194
<b>UniProt ID:</b>	<a href="#">Q8R550</a>
<b>Cytogenetics:</b>	X F4
<b>MW:</b>	69.3 kDa
<b>Gene Summary:</b>	Adapter protein involved in regulating diverse signal transduction pathways. Involved in the regulation of endocytosis and lysosomal degradation of ligand-induced receptor tyrosine kinases, including EGFR and MET/hepatocyte growth factor receptor, through an association with CBL and endophilins. The association with CBL, and thus the receptor internalization, may be inhibited by an interaction with PDCD6IP and/or SPRY2. Involved in regulation of ligand-dependent endocytosis of the IgE receptor. Attenuates phosphatidylinositol 3-kinase activity by interaction with its regulatory subunit. May be involved in regulation of cell adhesion; promotes the interaction between TTK2B and PDCD6IP. May be involved in the regulation of cellular stress response via the MAPK pathways through its interaction with MAP3K4. Is involved in modulation of tumor necrosis factor mediated apoptosis. Plays a role in the regulation of cell morphology and cytoskeletal organization. Required in the control of cell shape and migration (By similarity).[UniProtKB/Swiss-Prot Function]