

## Product datasheet for **MR226336L1V**

### Hnrnpab (NM\_010448) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hnrnpab (NM_010448) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hnrnpab
Synonyms:	3010025C11Rik; CBF-A; Cgbfa; Hnrpab
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_010448
ORF Size:	855 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226336).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_010448.3</a> , <a href="#">NP_034578.1</a>
RefSeq Size:	2404 bp
RefSeq ORF:	858 bp
Locus ID:	15384
UniProt ID:	<a href="#">Q99020</a>
Cytogenetics:	11 B1.3



[View online »](#)

**Gene Summary:**

This gene encodes a protein with consensus RNA binding domains present in a number of other RNA binding proteins and a glycine-rich C-terminus. This gene overlaps in a tail-to-tail orientation the gene encoding alanine-glyoxylate aminotransferase 2-like 2. Some of the exons of this gene are interspersed with exons of alanine-glyoxylate aminotransferase 2-like 2. Two alternatively spliced transcript variants that encode distinct proteins have been described for this gene. [provided by RefSeq, Jul 2008]