

Product datasheet for **MR226346**

Slco1a1 (NM_013797) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slco1a1 (NM_013797) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Slco1a1
Synonyms:	A530084B21; OATP-1; Oatp1; Oatp1a1; Slc21a1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR226346 representing NM_013797
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGAAGAAACAGAGAAAAGTTGCAACACAAGAAGGCAGATTCTTCTCCAAGATGAAGGTGTTTCTGA
 TGCATTAAACATGTGCATACCTAGCCAAATCACTATCAGGAGTTTATATGAATTCATGCTCACACAAAT
 AGAGAGACAATTTGGTATCCCTACATCTGTAGTTGGATTATCACTGGGAGCTTTGAGATAGGAAACCTT
 TTGTTGATTGTATTCGTGAGTTATTTTGAAGGAAACTGCACAGACCTATCATTATTGGTGTGGATGTG
 TGGTTATGGGCCTGGGGTGTTCCTTAATGGCATCACCTCATTTCCTCATGGGCAGATACAAATATGAAAC
 AACAAATACACCTACAAGCAACTGTCCTCAAACAGCTTCTTGTGTATAGAAAACAGAACACAGACCTTA
 AAGCCAACGCAAGATCCAACAGAGTGTGTAAAAGAAATAAAATCATTAAATGTGGATATATGTTACTAATAG
 GAAATACTATGCGTGAATTTGGTAAACTCCCATAAATGCCCTGGGTATTTCTTACATCGAAGACTTTGC
 CAAATCAGAAAACCTCCTTTATATATTGGGATTTTAGAAATGGGGAAGATTGTTGGCCCGATTATTGGA
 CTTTTGTTGGGATCTTTCTTTGCACGTGTTTATGTAGACATTGGGTCTGTGAATACAGATGACCTGACCA
 TAACTCCCCTGATACACGCTGGGTGCGTCTTGGTGGATTGGCTTTTTGGTCTGTGCAGGAGTGAATAT
 CCTGACCAGCATCCCCTTTTCTTTTCCAAAGACACTGCCAAAGAAAGAACTACAGGATAATGTGGAT
 GTGACTAAATATGAGAAAGTTGAAAACATAGAGAAAGAGCCAAAAGGAAAACCTTGGAACTCACTAAAG
 ATTTCTTGCCATTCATGAAGAGCCTCTGCTGCAATCCTATTTATATGCTGTTTCACTTACGAGTGTGCT
 CCAGATAAATGGATTTGCCAGTACATTTACCTTCTTGCCTTCCAGTATGCCTTGGATATTTAATCAGTG
 GCTTTATTATGAAGAAGTTAAGATCACTGTCAAGAAAGCTGCATACATAGCATTGGCCCTGCTTATC
 TGAATATTTTATTTTCTCTGTAACATTTGTTGACCTGTGATAATTTCCAGTTGCTGGATTAACCACC
 TCTTATAAAGGAGTTCAGCACCTCTATATGGGGAGAAAATGTCCTTGCTGACTGCAACACAAGGTGCA
 GCTGCTTAACGGACACATGGGATCCAGTGTGTGGGACAATGGCCTAGCTTACATGTCAGCCTGCCTCGC
 AGGCTGTGAGAAGTCTGTTGGAAGTGAACCAACATGGTGTTCAAAATGCAGCTGCATTGGGTGATCA
 GGAAACTCATCTGCAGTCTGGGCTGTGTAAGAAAGGCCCTGAGTGTGACAACAAGCTGCAGTACTTTT
 TAATCAAGTCAGTATTTCCAGTTTCATCTTCTCACTCGCAGCCATCCCTGGATACATGTTTCTTCTGAG
 GTGTGCAAACTGAAGAGAAGTCAATTGGAGTTGGATTACATGCATTTTTCATAAGACTATTAGCTGGC
 ATTCGGCACCTGTTTACTTTGGCGCTTTGATAGACAGAACCTGTTTACATTGGGGAAGCTGAAATGTG
 GTCAACCAGGGGCATGCAGGATGTATGATATAAACAGATTACAGGCACATTTACCTGGGGTTGCCTGCAGC
 AGTGAGAGGATCAAGCTTTCTCCCGCAGTCTTCAATCTAATCCTTATGAGGAAATTCATTTCCCTGGC
 GACATCCACTCTCCAGACACTGAACTTGACAGAGATGAAGCTCACAGAGAAGGAAAGCGAGTGCACAGATG
 GTGCAGAAGTCTAAGGTCGAGAACGATGGAGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGT

ACGCGTACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR226346 representing NM_013797
 Red=Cloning site Green=Tags(s)


MEETEKKVATQEGRRFFSKMKVFLMSLTCAYLAKSLSGVYMNSMLTQIERQFGIPTSVVGFITGSFEIGNL
 LLIVFVSYFGRKLHRPIIIGVGCVVMGLGCFLMASPHFLMGRYKYETTISPTSNLSSNSFLCIENRTQL
 KPTQDPTECVKEIKSLMWIYVLIGNTMRGIGETPIMPLGISYIEDFAKSENSPLYIGILEMGKIVGPIIG
 LLLGSFFARVYVDIGSVNTDDLITPTDTRWVGAWWIGFLVCAGVNILTSIPFFFFPKTLPKKELQDNVD
 VTKYEKVEKHRERAKKENLGITKDFLPFMKSLCCNPIYMLFSLTSVLQINGFASTFTFLPKYLEQQYGKS
 TSEAVFLIGVYSLPPVCLGYLISGFIMKKFKITVKAAYIAFGLSLSEYFIFLCNYLLTCDNFPVAGLTT
 SYKGVQHPPLYGEKNVLADCNTRCCLTDTWDPVCGDNGLAYMSACLAGCEKSVGTGTNMVFQNCSCIGSS
 GNSSAVLGLCKKGPECDNKLQYFLIKSVFSSFIFSLAAIPGYMVLRLCVKSEEKISIGVGLHAFFIRLLAG
 IPAPVYFGALIDRTCLHWGTLKCGQPGACRMYDINRFRHIYLGLPAAVRGSSFLPAVFILILMRKFHFPG
 DIHSPDTELAEMKLTEKESECTDVCSPKVENDGELKTKL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:




ACCN: NM_013797

ORF Size: 2010 bp

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_013797.5 , NP_038825.1
RefSeq Size:	3967 bp
RefSeq ORF:	2013 bp
Locus ID:	28248
UniProt ID:	Q9QXZ6
Cytogenetics:	6 G2
MW:	74.8 kDa
Gene Summary:	<p>Mediates the Na(+)-independent transport of organic anions such as taurocholate, prostaglandin E2 (PGE2), dehydroepiandrosterone sulfate (DHEAS), 17-beta-glucuronosyl estradiol, estrone-3-sulfate, sulfobromophthalein (BSP), ouabain and gadoxetate.</p> <p>[UniProtKB/Swiss-Prot Function]</p>

Product images:



Circular map for MR226346