

## Product datasheet for **MR226372L4V**

### Tial1 (NM\_009383) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tial1 (NM_009383) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tial1
Synonyms:	5330433G13Rik; AL033329; mTIAR; TIAR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_009383
ORF Size:	1176 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226372).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009383.2</a> , <a href="#">NP_033409.1</a>
RefSeq Size:	4292 bp
RefSeq ORF:	1179 bp
Locus ID:	21843
UniProt ID:	<a href="#">P70318</a>
Cytogenetics:	7 F3
Gene Summary:	RNA-binding protein. Possesses nucleolytic activity against cytotoxic lymphocyte target cells. May be involved in apoptosis (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)