

Product datasheet for **MR226373L4V**

Hoxc13 (NM_010464) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hoxc13 (NM_010464) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hoxc13
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010464
ORF Size:	984 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226373).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010464.2 , NP_034594.1
RefSeq Size:	2433 bp
RefSeq ORF:	987 bp
Locus ID:	15422
UniProt ID:	P50207
Cytogenetics:	15 57.99 cM
Gene Summary:	Transcription factor which plays a role in hair follicle differentiation. Regulates FOXQ1 expression and that of other hair-specific genes.[UniProtKB/Swiss-Prot Function]



[View online »](#)