

Product datasheet for **MR226397A1V**

Mouse Srsf9 (NM_025573) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Srsf9 (NM_025573) AAV Particle
Tag: Myc-DDK
Symbol: Srsf9
Synonyms: 25kDa; 2610029M16Rik; Sf; Sfrs9; SRp; SRp30c
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR226397 representing NM_025573
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGTCGTCGGGCTGGGCGGACGAGCGCGGCGGCGAGGGCGACGGGCGCATCTACGTGGGCAACCTTCGTC
 CCGACGTGCGCGAGAAGGACCTCGAGGACTTGTCTACAAGTACGGCCGCATCCGCGAGATCGAGCTCAA
 GAACCGGCACGGCCTCGTGCCCTTCGCCTTCGTGCGCTTCGAGGACCCGCGAGATGCTGAGGATGCGATC
 TATGGAAGAAACGGTTACGATTATGGCCAGTGTGACTCCGTGTGGAGTTCGCCAGGACTTACGGAGGTC
 GGGGTGGGTGGCCCGTGGTGCAAGGAACGGGCCTCTACAAGACGGTCAGATTTCGAGTTCCTGTTTC
 AGGACTTCCTCCATCAGGCAGCTGGCAGGACCTGAAAGATCACATGCGAGAAGCTGGGGATGCTGTATT
 GCAGACGTACAGAAGGACGGAATGGGGATGGTTGAATATTTGAGAAAAGAGGACATGGAATATGCTCTGC
 GTAAACTGGATGACACCAAATTCGCTCTCACGAGGTGAGACTTCTACATCCGAGTGTATCCTGAGAG
 AAGCACCAGCTATGGCTACTCACGGTCGCGGTCTGGGTCCAGGGGCCGCGACTCGCCATACCAAAGCCGG
 GGCTCGCCCACTACTTCTCTCTTTCAGGCCCTAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR226397 representing NM_025573
 Red=Cloning site Green=Tags(s)

MSSGWADERGGEGDGRIYVGNLPSDVREKDLDFYKYGRIREIELKNRHGLVPFAFVRFEDPRDAEDAI
 YGRNGYDYGQCLRVFPRTYGGRGGWPRGARNGPPTRRSDFRVLVSGLPSPGSWQDLKDHMREAGDVCI
 ADVQKDGMMGMEYLRKEDMEYALRKLDITKFRSHEGETSYIRVYPERSTSYGYSRSGSRGRDSPYQSR
 GSPHYFSPFRPY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_025573
ORF Size:	666 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_025573.3 , NP_079849.1
RefSeq Size:	1162 bp
RefSeq ORF:	669 bp
Locus ID:	108014
UniProt ID:	Q9D0B0
Cytogenetics:	5 F
MW:	26.1 kDa