

Product datasheet for **MR226399**

Srsf11 (NM_001093752) Mouse Tagged ORF Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | Srsf11 (NM_001093752) Mouse Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | Srsf11 |
| Synonyms: | 0610009J05Rik; 2610019N13Rik; BF642805; Sfr; Sfrs11 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)

ORF Nucleotide Sequence:

>MR226399 representing NM_001093752
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGTCGTACCGCCAGCTCCCGGGACCGCAGGACCGACGAGCCCTCCGTTGCCCTACCAACCTCGCCC
 TCAGGAGAAACTTCTTTTCGGGAGGCAGCCCGCCCATGAGCAGCACGGCAGTGGTCCCCAGCGCGCCAGG
 CCCGGGCCAGGCCCCAGTGGCGGGCCAGGCGGGCACCAGGTGATCCAGGTACCAACGTGTCTCCG
 AGCGCCAGCTCAGAGCAGATGCGTACGCTCTTCGGCTTCTGGCAAGATCGACGAGCTGCGCTCTTCC
 CGCCGGATGACTCACCTTGGCAGTCTCTCTCGTGTCTGCTTTGTTAAGTTCCATGATCCAGACTCAGC
 AGTTGTGGCACAGCATCTGACAAACACTGTGTTGTTGACAGAGCTTTGATAGTCGTACCATATGCTGAA
 GGAGTTATTCCTGATGAGACTAAGGCTTGTCTCTATTGGCACCAGCTAATGCAGTGGCAGGTCTTCTGC
 CTGGTGGTGGACTCTGCCTACCCCAACCCACTTACCCAGATTGGCGCTTCCCCTGGCTGCGTTGGG
 AGCTCCAGCTCTTGATCCTGCCCTTGCTGCTCTGGGCTTCCAGGAACAACTTGAAGTCTCAGTCCCTT
 GCTGCAGACCAGCTATTGAACTTATGAGTACTGTTGATCCCAAATTGAATCATGTAGCAGCTGGTCTTG
 TTTACCAAGTCTGAAGTCAGACACCTCTAGTAAAGAAATAGAGGAAGCCATGAAGAGAGTGCCTGAAGC
 CCAGTCCCTGATTTCTGCCGCTATTGAGCCAGACAAGAAGGAAGAAAAGAGAAGACATTCGAGATCAAGA
 TCTCGCTCGAGGAGGAGGAGGCCCCCTCATCTCCAGACACAGGCGGTACGAAGCCGTCAGGAGGC
 GGTCACTCCAATCCAGGAGTAGAAGGAGATCCAAAAGTCCAAGACGGAGACGATCTCATTCCCGAGA
 GAGGGGTAGAAGGTCAAGGAGCACATCCAAAGCCAGAGACAAAAGAAAGATAAAGAAAAGAAACGT
 TCCAAAACCTCCACAAAAGCTACAGCACAGCCAGACGCTCTCGAGTGAAGCAGAGAGAGAAGACGAC
 GGCGAAGTAGGAGTGGTACAAGATCCCTAAAAAGCCAGGTCTCCTAAAAGAAAGCTGTCTCGTCCACC
 ATCCCAAGGAGACATAAAAAGAGAAGAAGAAAGATAAGGACAAAAGAAAGCAGAGATGAAAGGGAA
 CGATCAACAAGTAAAAGAGAGGAGCAAGACAAAGGAGAAGGAGCGAGAAGGAAGTCTGAGAGTGACA
 AAGACGTCAAACAGGTGACTCGGGATTATGATGAAGAGGAGCAAGGCTATGACAGTGAGAAGAAAAGAA
 AGAGGAGAAGAGACCAACAGAAGCAGTGTCCCCTAAAACAAGGAGTGTCTGTGGAGAAGGGGTGGGC
 GACCTTCGAGAGTCAAAGTGAACGGGGATGATCACCACGAAGAAGACATGGATATGAGTGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR226399 representing NM_001093752
 Red=Cloning site Green=Tags(s)

MSSPPAPGTAGPTSPSVALTNLALRRNFFREAAAMSSTAVVPSAPGPGPGSPGGPGGTEVIQVTNVSP
 SASSEQMRTLFGFLGKIDELRLFPPDPSPLPVSSRVCFVKFHDPDSAVVAQHLTNTVFVDRALIVVPYAE
 GVIPDETKALSL LAPANAVAGLLPGGGLLPTPNPLTQIGAVPLAALGAPALDPAL AALGLPGTNLNSQSL
 AADQLLKMSTVDPKLNHVAAGLVSPSLKSDTSKEIEEAMKRVREAQSLISAAIEPDKKEEKRRHSR
 SRSRRRRTPSSSRHRSRHSRSHSKSRRRRSKSPRRRSHSRERGRSRSTSKARDKKKEDKEKKR
 SKTPPKSYSTARSRASRERRRRSRGTRSPKPKPRSPKRLSRSPSPRRHKKEKKDKDKERSRDERE
 RSTSKKRSKDKERERKSESDKVKQVTRDYDEEEQGYDSEKEKKEKRPTEAVSPKTKESVEKGVG
 DLRESKVNDDHHEEDMDMSD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001093752

ORF Size: 1533 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001093752.3](#)

RefSeq Size: 3086 bp

RefSeq ORF: 1536 bp

Locus ID: 69207


Cytogenetics: 3 H4

MW: 57.1 kDa

Gene Summary:

The protein encoded by this gene is a member of the serine/arginine (SR)-rich family of pre-mRNA splicing factors, which constitute part of the spliceosome. Each of these factors contains an RNA recognition motif (RRM) for binding RNA and an RS domain for binding other proteins. The RS domain is rich in serine and arginine residues and facilitates interaction between different SR splicing factors. In addition to being critical for mRNA splicing, the SR proteins have also been shown to be involved in mRNA export from the nucleus and in translation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2014]

Product images:



Circular map for MR226399