

## Product datasheet for **MR226453L4V**

### Usp27x (NM\_019461) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Usp27x (NM_019461) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Usp27x
Synonyms:	Sfc11; Usp27
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_019461
ORF Size:	1314 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226453).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_019461.4</a> , <a href="#">NP_062334.2</a>
RefSeq Size:	3240 bp
RefSeq ORF:	1317 bp
Locus ID:	54651
UniProt ID:	<a href="#">Q8CEG8</a>
Cytogenetics:	X 3.31 cM
Gene Summary:	Deubiquitinase that can reduce the levels of BCL2L11/BIM ubiquitination and stabilize BCL2L11 in response to the RAF-MAPK-degradation signal. By acting on BCL2L11 levels, may counteract the anti-apoptotic effects of MAPK activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)