

Product datasheet for **MR226454**

Eif4a1 (NM_001159375) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Eif4a1 (NM_001159375) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Eif4a1
Synonyms: BM-010; Ddx2a; Eif4
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR226454 representing NM_001159375
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTCTGCGAGTCAGGATTCTCGATCCAGAGACAATGGCCCCGACGGGATGGAGCCGGAAGGCGTCATCG
 AGAGTAACTGGAACGAGATTGTGGATAGCTTTGATGACATGAATCTCTCAGAGTCCCTCCTCGTGGTAT
 TTATGCCTATGGTTTTGAGAAGCCCTCTGCCATCCAGCAGCGAGCTATTCTTCTGTATCAAGGGTTAT
 GATGTGATTGCTCAAGCCAGTCTGGGACTGGGAAAACAGCTACATTTGCCATATCAATTCTGCAGCAGA
 TTGAATTAGATCTAAAGCCACTCAGGCTTTGGTTCTGGCACCCACACGTGAATTGGCTCAGCAGATACA
 AAAGGTGGTTATGGCATTAGGAGACTACATGGGTGCCTCTTGTCATGCCTGCATTGGGGGCACCAATGTG
 CGTGCTGAGGTGCAGAAGCTGCAGATGGAAGCTCCCCATATCATCGTGGGTACCCCTGGCCGGGTGTTTG
 ACATGCTTAACCGGAGATACCTGTCCCCAAATACATCAAGATGTTCTGACTGGATGAAGCAGATGAAAT
 GTTAAGCCGAGGGTTCAAGGATCAGATCTATGACATATTCCAGAAGCTCAACAGCAACACACAGGTAGTT
 TTGTTGCTGCTACAATGCCTTCTGATGTCTTGAGGTGACCAAGAAATTTATGAGAGCCCTATTTCGGA
 TTCTTTGCAAGAAGGAAGAATTGACCCTGGAGGGTATCCGCCAATTCTACATCAATGTGGAACGAGAGGA
 GTGGAAGCTTGACACATTGTGTGACTTGTATGAGACGCTGACCATCACCCAGGCAGTCATCTTTATCAAC
 ACCAGAAGGAAGGTGGACTGGCTCACCGAGAAGATGCATGCCCGAGATTTCACTGTTTCTGCCATGCACG
 GAGATATGGACAAAAGGAACGAGATGTGATCATGAGGGAGTTCCGGTCTGGCTCTAGCAGAGTATTAAT
 TACCACTGACCTGTTGGCCAGAGGCATTGATGTGCAGCAGGTCTCCTTAGTCATCAACTATGACCTTCCC
 ACCAACAGGGAAAACACTACATCCACAGGTGGTTCGGTTTGGTTCG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >MR226454 representing NM_001159375
 Red=Cloning site Green=Tags(s)

MSASQDSRSRDNGPDGMEPEGVIESNWNEIVDSFDDMNLSESLLRGIYAYGF EKPSAIQQRILPCIKGY
 DVIAQAQSGTGKTATFAISILQQIELDLKATQALVLAPTRELAQQIQKVMALGDYMGASCHACIGGTNV
 RAEVQKLQMEAPHIIVGTPGRVFDMLNRRYLSPKYIKMFVLDEADEMLSRGFKDQIYDIFQKLSNTQVV
 LLSATMPSDVLEVTKKFMRDPIRILVKKEELTEGIRQFYINVEREEWKLDLTLCDLYETLTITQAVIFIN
 TRRKVDWLTEKMHARDFTVSAMHGMDQKQERDVIMREFRSGSSRVLITDLLARGIDVQVSLVINYDLP
 TNRENYIHRWSVWS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001159375

ORF Size:	1092 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001159375.1 , NP_001152847.1
RefSeq Size:	1886 bp
RefSeq ORF:	1095 bp
Locus ID:	13681
Cytogenetics:	11 42.86 cM
MW:	41.9 kDa
Gene Summary:	ATP-dependent RNA helicase which is a subunit of the eIF4F complex involved in cap recognition and is required for mRNA binding to ribosome. In the current model of translation initiation, eIF4A unwinds RNA secondary structures in the 5'-UTR of mRNAs which is necessary to allow efficient binding of the small ribosomal subunit, and subsequent scanning for the initiator codon.[UniProtKB/Swiss-Prot Function]