

Product datasheet for **MR226464L3V**

Carns1 (NM_134148) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Carns1 (NM_134148) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Carns1
Synonyms:	AI790298; Atpgd1; mKIAA1394
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_134148
ORF Size:	2841 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226464).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_134148.2 , NP_598909.2
RefSeq Size:	3968 bp
RefSeq ORF:	2844 bp
Locus ID:	107239
Cytogenetics:	19 A
Gene Summary:	Catalyzes the synthesis of carnosine and homocarnosine. Carnosine is synthesized more efficiently than homocarnosine.[UniProtKB/Swiss-Prot Function]



[View online »](#)