

Product datasheet for **MR226466L3V**

Eif4a2 (NM_001123038) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Eif4a2 (NM_001123038) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Eif4a2
Synonyms:	4833432N07Rik; BM-010; Ddx2b; eIF-4A-II; Eif4; eIF4A-II
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001123038
ORF Size:	1086 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226466).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001123038.1 , NP_001116510.1
RefSeq Size:	3174 bp
RefSeq ORF:	1089 bp
Locus ID:	13682
UniProt ID:	P10630
Cytogenetics:	16 13.93 cM



[View online »](#)

Gene Summary:

ATP-dependent RNA helicase which is a subunit of the eIF4F complex involved in cap recognition and is required for mRNA binding to ribosome. In the current model of translation initiation, eIF4A unwinds RNA secondary structures in the 5'-UTR of mRNAs which is necessary to allow efficient binding of the small ribosomal subunit, and subsequent scanning for the initiator codon.[UniProtKB/Swiss-Prot Function]