

Product datasheet for **MR226471L3V**

Krt71 (NM_019956) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Krt71 (NM_019956) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Krt71
Synonyms:	AA589543; Ca; Cal4; Cu; Krt2-6g; mK6irs; mK6irs1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019956
ORF Size:	1575 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226471).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019956.1 , NP_064340.1
RefSeq Size:	2193 bp
RefSeq ORF:	1575 bp
Locus ID:	56735
UniProt ID:	Q9R0H5
Cytogenetics:	15 57.0 cM
Gene Summary:	Plays a central role in hair formation. Essential component of keratin intermediate filaments in the inner root sheath (IRS) of the hair follicle.[UniProtKB/Swiss-Prot Function]



[View online »](#)