

## Product datasheet for **MR226485L3V**

### Usp37 (NM\_176972) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Usp37 (NM_176972) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Usp37
Synonyms:	4932415L06Rik; C330008N13Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_176972
ORF Size:	2940 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226485).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_176972.4</a> , <a href="#">NP_795946.2</a>
RefSeq Size:	7196 bp
RefSeq ORF:	2940 bp
Locus ID:	319651
UniProt ID:	<a href="#">Q8C0R0</a>
Cytogenetics:	1 C3



[View online »](#)

**Gene Summary:**

Deubiquitinase that antagonizes the anaphase-promoting complex (APC/C) during G1/S transition by mediating deubiquitination of cyclin-A (CCNA1 and CCNA2), thereby promoting S phase entry. Specifically mediates deubiquitination of 'Lys-11'-linked polyubiquitin chains, a specific ubiquitin-linkage type mediated by the APC/C complex. Also mediates deubiquitination of 'Lys-48'-linked polyubiquitin chains in vitro. Phosphorylation at Ser-628 during G1/S phase maximizes the deubiquitinase activity, leading to prevent degradation of cyclin-A (CCNA1 and CCNA2). Plays an important role in the regulation of DNA replication by stabilizing the licensing factor CDT1.[UniProtKB/Swiss-Prot Function]