

Product datasheet for **MR226524L3V**

Tnfrsf13c (NM_028075) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tnfrsf13c (NM_028075) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tnfrsf13c
Synonyms:	2010006P15Rik; BAFF-R; Baffr; Bcmd; Bcmd-1; Bcmd1; Lvis22
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_028075
ORF Size:	528 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226524).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_028075.2
RefSeq Size:	1906 bp
RefSeq ORF:	528 bp
Locus ID:	72049
UniProt ID:	Q9D8D0
Cytogenetics:	15 E1
Gene Summary:	B-cell receptor specific for TNFSF13B/TALL1/BAFF/BlyS. Promotes the survival of mature B-cells and the B-cell response.[UniProtKB/Swiss-Prot Function]



[View online »](#)