

Product datasheet for **MR226525**

Bcl11b (NM_001079883) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Bcl11b (NM_001079883) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Bcl11b
Synonyms:	9130430L19Rik; AI604821; B630002E05Rik; BCL-11B; Ctip2; Rit1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR226525 representing NM_001079883
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATG**TCCCGCCGCAAA**CAGGGCA**ACC**CGCAGCACTTGTCCAGAGGGA**ACT**CATCACGCCAGAGGCTGACC
 ATGTGGAGGCTACCATCCTCGAGGAAGACGAGGGTCTGGAGATAGAGGAGCCTAGCAGCCTGGGGCTGAT
 GGTGGGAGGCCCGACCTGATCTACTCACCTGTGGCCAGTGTGAGATGAACTTCCCGCTGGGGACATC
 CTGTTTTTATAGAGCACAAGAAGAAACAGTGTGGAGGCTGGGCCCTGCTACGACAAGGTCTGGACA
 AGAGCAGTCCACCTCCCTCCTCTCGCTCTGAGCTCAGGAGAGTATCTGAGCCAGTGGAGATCGGGATCCA
 GGT**CACCC**CTGATGAAGATGACCACCTACTGTCACCCACGAAAGGCATCTGTCCAAGCAGGAGA**ACATT**
 GCAGGGCCGTGCAGGCC**TGCC**AGCTGCCATCGATGGCCCCATAGCTGCCTCCTTCCACCC**TCCCA**
 CCTCCGTGATTACTT**CACCT**CTGCGTGGCCTGGGGCTTCTCCACCC**TGTTCCCGCTGCCTT**GCTGTGG
 TGCACGCC**CCATCTCGGGC**ACGGGACTCAGGGT**GAGGGT**CAGATGGAGGCTCC**TTTGGATGCC**AGTGT
 GAGTTGT**CAGGTA**AAAGATGAGCCTTCCAGCTACATTTGCACAACATGCAAGCAGCC**TTCAAC**AGCGCCT
 GGTTCTGCTGCAGCACGCACAGAACACACATGGCTTCCGAATCTACCTGGAGCCTGGGGCGCCAGCAC
 CTCGCTCACGCCAGGCTCACCATCCCGCCACC**GCTCGGG**CGGAGACCGTGGCGCAGTCC**CCACT**CATG
 AATTTCTGGGGACAGCAATCCTTTCAACCTGCTGCGCATGACGGG**CCCATCCTGCGG**ACCAC**CTG**
 GCTTCGGT**GAGGG**CCGTTGCCAGGTACGCCACCGCTTTCAGCCACC**GCCACGCC**ATCACTTGGACCC
 ACACCGCCTCAGTGCAGAGGAGATGGGGCTCGTGGCCAGCACCC**AGTGCC**TT**CGAC**CGAGT**CATGCGC**
 CTGAACCC**ATGGCC**ATAGACTCTCCTGCCATGGACTTCTCCCGCGGCTGCAGAG**ACTGGCC**CAACA
 GCTCCAGCCCGCCCGCTGTCCCAGGCGTGGCAACC**ATGCA**CCGGT**GCTGA**ACC**TTTCC**AGCC
 CAGTCCCAAGTCCCGTTCTCAGCACGCCACC**GCTGCC**ACC**ATGCCTGCGGG**ACACC**GCCACCG**CAG
 CCGCTGCCAAGAGCAAGTCTGTGAGTCTGCGGCAAGAC**TTCAAGT**CCAGAGCAATCTCATCGTGC
 ACCGGCGCAGCCACACGGGCGAGAAGCC**CTACA**AGTGCCAGTGTGCGACCATGCGTGTGCGAGGCGAG
 CAAGCTCAAGCGCCACATGAAGACGCACATGCACAAGGCGGGCTCTTGCTGGCCGCTCAGACGACGGG
 CTCTCAGTGCAGCTCCCTGAGCCGGGACCAGCGAGCTGCCAGGTGACCTGAAAGCGGCCGATGGCG
 ACTTCCGCCACCATGAGAGCGACCCATCTTGGCC**CCGAGCCTGAGG**ACGACGAGGACGAGGAGGAGGA
 AGAAGAGGAGCTGCTGCTGGAGAACGAGAGCCGGCCTGAGTCGAGCTT**CAG**CATGGACTCGGAGCTGGGC
 CGTGGCCGCGAGAACGGAGGTGGCGTGCACCGGGGGTGGCGGGCGCAGGGGCTGCAGCTGCGGCTCTGG
 CGGATGAGAAGGCTCTGGCCCTGGGCAAGGTGATGGAGGACGCAGGGCTGGGCGCACTGCCG**CAGTATGG**
 GGAGAAGCGGGGCGCCTTCTGAAGCGTGCAGGCGACCGGGTATGCCGGAGCTGTTGGCTGTGGGGAC
 GCGGGTGCACCGGGTGCAGTGAACGGGCGCGGGGGCC**TTCGCG**CCAGGGCGAGAGCC**TTTCC**AGCTC
 TCTTCCACGCAAGCCAGCACCGCTGCCAGCC**TGGGCTCGGTGGT**CCCGCGCTGCACGCGGCCAAGCG
 CATCAAGGTGGAGAAAGACCTGGAGCTGCCACCTGCCGCCCTCATCCCATCTGAGAACGTGTACTCGCAG
 TGGCTCGTGGGCTACGCAGCATCGCGCCACTT**CATGA**AGGACCCATTCTGGGCTT**CACGG**ATGCGCGCC
 AGTGCCTTT**CGCC**ACATCGTCGGAACATTCTCTGAGAACGGCAGCCTGCGCTTCTCAACGCCACCCGG
 GGACTGCTGGACGGCGGGCTGTCCGGGCGAGTGGCACGGCGAGCGGGGCGAGCACACCTCACCTGGGT
 GGTCCGGTCTGGGAGGCCGAGCTCCAAGGAGGGCCGCGCAGCGACACATGTGAGTACTGCGGCAAGG
 TCTTCAAGAACTGTAGCAACCTGACGGTGCACCGGAGGAGCCACACCGGCGAGCGGCCTTACAAGTGCGA
 GCTGTGCAACTACGCGTGCAGCGCAGAGCAGCAAGCTCACGCGCCACATGAAGACGCACGGGCGAGATCGGC
 AAGGAGGTGTACCCTGCGACATCTGCCAGATGCCCTT**CAG**CGCTACAGCACCC**TGGAG**AAACACATGA
 AAAAGTGGCACGGTGAACACTT**GCTGACT**AATGATGTCAAATCGAGCAGGCTGAGAGGAGC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR226525 representing NM_001079883
 Red=Cloning site Green=Tags(s)


MSRRKQGNPQHLSQRELITPEADHVEATILEEDEGLEIEEPSSLGLMVGPPDLLTCGQCQMNFLGDI
 LVFIEHKKKQCGLGPCYDKVLDKSSPPSSRSELRRVSEPVVEIGIQVTPDEDDHLLSPTKIGICPKQENI
 AGPCRPAQLPSMAPIAASSSHPPTSIVITSPRLALGVLPPCFPLPCCGARPISGDGTQEGEQMEAPFGCQC
 ELSGKDEPSSYICTTCKQPFNSAWFLLQHAQNTGHFRIYLEPGPASTSLTPRLTIPPLGPETVAQSPLM
 NFLGDSNPFNLLRMTGPIILRDHPGFEGRLPGTPPLFSPPPRHLLDPHRLSAEEMGLVAQHPSAFDRVMR
 LNPMAIDSPAMDFSRRLRELAGNSSTPPVSPGRGNPMHRLLNPFQSPKSPFLSTPPLPMPAGTPPPQ
 PPAKSKSCEFCGKTFKFQSNLIVHRSHTGEKPYKCQLCDHACSQASKLKRHMKTHMHKAGSLAGRSDDG
 LSAASSPEPGTSELPGLKAADGDFRHHESDPSLGPEPEDEDEEEEEELLENESRPESFSMDSELG
 RGRENGGGVPPGVAGAGAAAAALADEKALALGKVMEDAGLGALPQYGEKRGAFKRAAGTGDAGAVGCGD
 AGAPGAVNRRGGAFAPGAEPFALFPRKPAPLPSPLGGPALHAAKRIKVEKDLELPPAALIPSENVYSQ
 WLVGAAASRHFMKDPFLGFTDARQSPFATSSEHSENGSLRFSTPPGDLLDGGLSGRSGTASGGSTPHLG
 GPGGPRPSSKEGRRSDTCEYCGKVFKNCSNLTVHRSHTGERPYKCELCNYACAQSSKLRHMKTHGQIG
 KEVYRCDICQMPFSVYSTLEKHMKKWHGEHLLTNDVKIEQAERS

SGP TRTRRL EQKLI SEEDLA ANDILDYK DDDDKV

Restriction Sites:

SgfI-RsrII


Cloning Scheme:



* The last codon before the Stop codon of the ORF

ACCN:	NM_001079883
ORF Size:	2652 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001079883.1 , NP_001073352.1
RefSeq Size:	7923 bp
RefSeq ORF:	2655 bp
Locus ID:	58208
UniProt ID:	Q99PV8
Cytogenetics:	12 59.1 cM
MW:	95 kDa
Gene Summary:	Key regulator of both differentiation and survival of T-lymphocytes during thymocyte development in mammals (PubMed:12717433). Essential in controlling the responsiveness of hematopoietic stem cells to chemotactic signals by modulating the expression of receptors CCR7 and CCR9, which direct the movement of progenitor cells from the bone marrow to the thymus (By similarity). Is a regulator of IL2 promoter and enhances IL2 expression in activated CD4(+) T-lymphocytes (PubMed:16809611). Tumor-suppressor protein involved in T-cell lymphomas. May function on the P53-signaling pathway. Repress transcription through direct, TFCOUP2-independent binding to a GC-rich response element.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR226525