

# Product datasheet for MR226525L3

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Bcl11b (NM\_001079883) Mouse Tagged Lenti ORF Clone

### **Product data:**

**Product Type:** Expression Plasmids

**Product Name:** Bcl11b (NM\_001079883) Mouse Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: Bcl11b

**Synonyms:** 9130430L19Rik; Al604821; B630002E05Rik; BCL-11B; Ctip2; Rit1

Mammalian Cell Puromycin

Selection:

**Vector:** pLenti-C-Myc-DDK-P2A-Puro (PS100092)


E. coli Selection: Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(MR226525).

Sequence:

Restriction Sites: Sgfl-Rsrll

**Cloning Scheme:** 






<sup>\*</sup> The last codon before the Stop codon of the ORF.





#### Plasmid Map:



**ACCN:** NM\_001079883

ORF Size: 2652 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001079883.1, NP 001073352.1</u>

RefSeq Size: 7923 bp RefSeq ORF: 2655 bp



## Bcl11b (NM\_001079883) Mouse Tagged Lenti ORF Clone - MR226525L3

 Locus ID:
 58208

 UniProt ID:
 Q99PV8

 Cytogenetics:
 12 59.1 cM

**Gene Summary:** Key regulator of both differentiation and survival of T-lymphocytes during thymocyte

development in mammals (PubMed:12717433). Essential in controlling the responsiveness of hematopoietic stem cells to chemotactic signals by modulating the expression of receptors CCR7 and CCR9, which direct the movement of progenitor cells from the bone marrow to the thymus (By similarity). Is a regulator of IL2 promoter and enhances IL2 expression in activated CD4(+) T-lymphocytes (PubMed:16809611). Tumor-suppressor protein involved in T-cell

lymphomas. May function on the P53-signaling pathway. Repress transcription through direct, TFCOUP2-independent binding to a GC-rich response element.[UniProtKB/Swiss-Prot

Function]