

Product datasheet for **MR226533L4V**

P2ry13 (NM_028808) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	P2ry13 (NM_028808) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	P2ry13
Synonyms:	2010001L06Rik; GPCR1; Gpr86; GPR94; P2Y13; P2yr13; SP174
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_028808
ORF Size:	1014 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226533).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_028808.3 , NP_083084.2
RefSeq Size:	2541 bp
RefSeq ORF:	1014 bp
Locus ID:	74191
UniProt ID:	Q9D8I2
Cytogenetics:	3 D
Gene Summary:	Receptor for ADP. Coupled to G(i)-proteins. May play a role in hematopoiesis and the immune system (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)