

Product datasheet for **MR226561L3V**

Meis2 (NM_010825) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Meis2 (NM_010825) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Meis2
Synonyms:	A430109D20Rik; Mei; Mrg; Mrg1; Str; Stra10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010825
ORF Size:	1410 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226561).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010825.3 , NP_034955.1
RefSeq Size:	4621 bp
RefSeq ORF:	1413 bp
Locus ID:	17536
UniProt ID:	P97367
Cytogenetics:	2 58.28 cM



[View online »](#)

Gene Summary:

This gene encodes a homeobox protein belonging to the TALE ('three amino acid loop extension') family of homeodomain-containing proteins. TALE homeobox proteins are highly conserved transcriptional regulators and several members have been shown to be essential contributors to developmental programs. In mice, a knock-out of this gene leads to lethality at embryonic day 14, accompanied with hemorrhaging. Embryos lacking this gene show defects in tissues derived from the neural crest, suggesting a critical role of this gene during cranial and cardiac neural crest cell development. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2016]