

## Product datasheet for **MR226562L3V**

### Sulf2 (NM\_028072) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Sulf2 (NM_028072) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sulf2
Synonyms:	2010004N24Rik; AU020235; mKIAA1247; MSulf-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_028072
ORF Size:	2625 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226562).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_028072.4</a> , <a href="#">NP_082348.2</a>
RefSeq Size:	3790 bp
RefSeq ORF:	2628 bp
Locus ID:	72043
UniProt ID:	<a href="#">Q8CFG0</a>
Cytogenetics:	2 H3
Gene Summary:	Exhibits arylsulfatase activity and highly specific endoglucosamine-6-sulfatase activity. It can remove sulfate from the C-6 position of glucosamine within specific subregions of intact heparin (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)