

## Product datasheet for **MR226587L3V**

### Ctps2 (NM\_001168571) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ctps2 (NM_001168571) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ctps2
Synonyms:	A830031M15Rik; A1326475; Ctpsh
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001168571
ORF Size:	1605 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226587).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001168571.1</a> , <a href="#">NP_001162043.1</a>
RefSeq Size:	5453 bp
RefSeq ORF:	1608 bp
Locus ID:	55936
UniProt ID:	<a href="#">P70303</a>
Cytogenetics:	X F4
Gene Summary:	Catalyzes the ATP-dependent amination of UTP to CTP with either L-glutamine or ammonia as the source of nitrogen. Constitutes the rate-limiting enzyme in the synthesis of cytosine nucleotides (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)