

Product datasheet for MR226592L4

Irf5 (NM 012057) Mouse Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Irf5 (NM_012057) Mouse Tagged Lenti ORF Clone

Tag: mGFP Symbol: Irf5

Synonyms: AW491843; mirf5

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(MR226592).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_012057

ORF Size: 1491 bp



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

Irf5 (NM_012057) Mouse Tagged Lenti ORF Clone - MR226592L4

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 012057.3</u>, <u>NP 036187.1</u>

 RefSeq Size:
 2081 bp

 RefSeq ORF:
 1494 bp

 Locus ID:
 27056

 UniProt ID:
 P56477

Cytogenetics: 6 12.36 cM


Gene Summary: Transcription factor involved in the induction of interferons IFNA and INFB and inflammatory

cytokines upon virus infection. Activated by TLR7 or TLR8 signaling (By similarity).


[UniProtKB/Swiss-Prot Function]



Product images:



Circular map for MR226592L4



Double digestion of MR226592L4 using Sgfl and Mlul