

Product datasheet for **MR226655A1V**

Mouse FasI (NM_010177) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse FasI (NM_010177) AAV Particle
Tag:	Myc-DDK
Symbol:	FasI
Synonyms:	APT1LG1; CD95-L; CD95L; CD178; Fas-L; Faslg; gld; Tnfsf6; Tnlg1a
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR226655 representing NM_010177 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCAGCAGCCCATGAATTACCCATGTCCCCAGATCTTCTGGGTAGACAGCAGTGCCACTTCATCTTGGG
CTCCTCCAGGGTCAGTTTTTCCCTGTCCATCTTGTGGGCTAGAGGGCCGGACCAAAGGAGACCGCCACC
TCCACCACCACCTGTGTCACTACTACCACCGCCATCACAACCACTCCCACTGCCGCCACTGACCCCTCTA
AAGAAGAAGGACCACAACACAAATCTGTGGCTACCGGTGGTATTTTTCATGGTTCTGGTGGCTCTGGTTG
GAATGGGATTAGGAATGTATCAGCTCTTCCACCTGCAGAAGGAAGTGGCAGAACTCCGTGAGTTACCAA
CCAAAGCCTTAAAGTATCATCTTTTGAAGCAAAATAGCCAACCCAGTACACCCCTCTGAAAAAAGAG
CCGAGGAGTGTGGCCATTTAACAGGGAACCCCACTCAAGGTCCATCCCTCTGGAATGGGAAGACACAT
ATGGAACCGCTCTGATCTCTGGAGTGAAGTATAAGAAAGGTGGCCTTGTGATCAACGAACTGGGTTGTA
CTTCGTGTATTCCAAAGTATACTTCCGGGGTCAGTCTTGCAACAACAGCCCTAAACCACAAGGTCTAT
ATGAGGAAGTCTAAGTATCCTGAGGATCTGGTGCTAATGGAGGAGAAGAGGTTGAAGTACTGCACTACTG
GACAGATATGGGCCACAGCAGTACCTGGGGCAGTATTCAATCTTACCAGTGTGACCATTTATATGT
CAACATATCTCAACTCTCTGATCAATTTGAGGAATCTAAGACCTTTTCGGCTTGATAAGCTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR226655 representing NM_010177 Red=Cloning site Green=Tags(s) MQQPMNYPCQIFWVDSSATSSWAPPGSVFPCPSCGPRGPDQRRPPPPPPVVSPLPPSQPLPLPPLTPL KKKDHNTNLWLPVVFVMVLVALVGMGLGMYQLFHLQKELAELEFNTQSLKVSSFEKQIANPSTPSEKKE PRSVAHLTGNPHSRISPLEWEDTYGTALISGVKYKGGVLINETGLYFVYSKVYFRGQSCNNQPLNHHKVY MRNSKYPEDLVLMEEKRLNYCTTGQIWAHSSYLGAVFNLTSADHLYVNIQSLSLINFEEKTFFGLYKL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_010177
ORF Size:	837 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_010177.4, NP_034307.1</u>
RefSeq Size:	1935 bp
RefSeq ORF:	840 bp
Locus ID:	14103
UniProt ID:	<u>P41047</u>
Cytogenetics:	1 69.95 cM
MW:	31.9 kDa