

## Product datasheet for **MR226657**

### Vldlr (NM\_013703) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Vldlr (NM_013703) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Vldlr
Synonyms:	AA408956; AI451093; AW047288
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>MR226657 representing NM\_013703  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**


ATGGGCACGTCCGCGCTGGGCGCTGTGGCTGCTGCTCGCGCTGTGCTGGGCGCCCCGGGACAGCGGCG  
 CCACTGCAAGCGGAAGAAGCCAAATGTGATAGCTCCAGTTTCAGTGACAAAATGGCCGCTGCATTAC  
 CCTGCTGTGAAATGTGATGGAGATGAAGACTGTGCGGATGGCAGCGACGAGAAGAACTGTGTAAGAAG  
 ACGTGTGCTGAGTCTGACTTCGTGTGCAAAAACGGCCAGTGTGTTCTAACAGATGGCAGTGTGACGGGG  
 ATCCTGATTGCGAAGACGGTTCTGATGAAAGCCCTGAACAGTGCCATATGAGAACATGCCGATAAATGA  
 AATCAGCTGTGGCGCCGTTCTACTCAGTGTATCCCCGTGCTCTGGAGATGCGATGGTGAATGATTGT  
 GACAATGGAGAAGATGAAGAAAATGTGGCAACATAACATGTAGTGCAGATGAGTTCCTGCTCCAGTG  
 GCCGCTGCGTCTCCAGAACTTTGTGTCAATGGCCAGGATGACTGTGACGATGGCAGTGTGAGCTGGA  
 CTGTGCTCCACCAACCTGCGGAGCCACGAGTTCAGTGCAGCACCTCTTCTGCATTCCCCTCAGCTGG  
 GTGTGTGATGATGACGCAGACTGTTAGACCAATCAGACGAGTCTTGGAGCAGTGTGGCGTCCAGCTG  
 TGATACATACCAATGTCTACCAAGTGAAGTCCAGTGTGGCTCTGGCGAGTGCATTACAAAAAATGGCG  
 GTGTGACGGAGACCCTGACTGCAAGGACGGCAGCGATGAGGTCAACTGCCCTTCTCGAACCTGCCGACCT  
 GACCAGTTGAATGTGAAGATGGTAGCTGTATCCACGGCAGCAGGCAATGCAATGGCATCCGAGACTGTG  
 TTGATGGCTCTGATGAAGTCAACTGCAAAAACGTCATCAGTGCCTGGGCCCTGGAAAGTTCAAGTGCAG  
 AAGCGGGAAATGCATAGACATGAGCAAAGTATGTGACCAGGAACAAGACTGCAGAGACTGGAGTGCAGG  
 CCCCTGAAGGAATGCCATATCAACGAATGCCTGGTCAATAATGGTGGCTGTTCCCATATGCAAAAGACC  
 TAGTTATAGGTTATGAGTGTGTTGTGACGCTGGGTTTGAAGTATGATAGGAAAACCTGTGGAGACT  
 TGATGAATGCCAAAACCGGGGATCTGCAGTCAAATTTGTATCAACTTAAAAGGCGGTTACAAGTGTGAA  
 TGTAGTCTGGCTATCAATGGATCTTGCCACTGGCGTGTGCAAGGCAGTAGGCAAAAGAGCCGAGTCTGA  
 TCTTCACTAATCGAAGAGACATCAGGAAGATTGGCCTAGAGAGAAAGGAATACATCCAATTGTAGAGCA  
 ACTAAGGAACACGGTGGCTCTCGATGCGGACATTGCAGCTCAGAAGCTGTTTTGGGCTGATCTCAGCCAG  
 AAGGCCATCTTCAAGTGCCTCAATTGATGACAAGTTGGTAGACATTTTAAAATGATCGACAATGTCTATA  
 ATCCTGCAGCCATTGCTGTTGATTGGGTGTACAAGACCCTACTGGACTGATGCGGCTTCTAAGACTAT  
 TTCAGTAGCTACCCTAGACGGAGCCAAGAGGAAGTTCTGTAAATCTGACTTGCAGAGCCTGCCTCC  
 ATAGCTGTGGATCCGTTGTCGGGCTTTGTTACTGGTCAAGTGGGCGAGCCAGCTAAAATAGAAAAAG  
 CAGGAATGAATGATTTGATAGACGTCCTCTGGTGACGGAGGACATCCAATGGCCTAATGGAATTACACT  
 CGACCTTGTCAAAAGCCGCTCTACTGGCTGGATTCCAAGTTCACATGTCTCTAGTGTGGACCTGAAT  
 GGTCAAGATCGTAGGATAGTGTCAAGTCTCTGGAGTTCCTAGCTCATCCTCTTGCACTCACCATATTTG  
 AGGATCGCGTCTACTGGATAGATGGAGAAAATGAAGCAGTGTACGGTGCCAATAAATTCCTGGGTCAGA  
 GCTGGCCACTCTAGTGAATAATCTCAATGATGCCAAGACATCATTGTCTACCATGAACCTCGTCCAGCCG  
 TCAGGTA AAAACTGGTGTGAAGACGATATGGAGAATGGAGGATGTGAATATCTCTGCCTGCCAGCACCAC  
 AGATCAATGACCACTCTCAAAAATACCTGTTCTGTCCCAATGGGTACAATCTCGAAGAAAATGGACG  
 AGAGTGTCAAAGTACTTCAACTCTGTGACTTACAGTGAGACAAAAGATATCAACACAACAGACATTCTA  
 CGAACTAGTGGACTGGTTCCTGGAGGGATCAATGTGACCACAGCAGTATCAGAAGTCAAGTGTCCCCAA  
 AAGGGACTTCAGCTGCCTGGGCCATCCTTCTCTTGTCTTAGTGTGCGAGCAGTAGGTGGCTACTT  
 GATGTGGAGGAATTGGCAACATAAAAACATGAAAAGCATGAACCTTGGACAATCCTGTGACTTGAAGACC  
 ACTGAAGAGGACCTGTCGATAGACATTGGTAGACACAGCGCTTCTGTAGGACACACATACCAGCAATAT  
 CAGTTGTAAGCACAGATGATGATCTGGCT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA



<b>ACCN:</b>	NM_013703
<b>ORF Size:</b>	2619 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_013703.2</a> , <a href="#">NP_038731.2</a>
<b>RefSeq Size:</b>	8327 bp
<b>RefSeq ORF:</b>	2622 bp
<b>Locus ID:</b>	22359
<b>UniProt ID:</b>	<a href="#">P98156</a>
<b>Cytogenetics:</b>	19 21.77 cM
<b>MW:</b>	96.8 kDa
<b>Gene Summary:</b>	Binds VLDL and transports it into cells by endocytosis. In order to be internalized, the receptor-ligand complexes must first cluster into clathrin-coated pits. Binding to Reelin induces tyrosine phosphorylation of Dab1 and modulation of Tau phosphorylation. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR226657