

## Product datasheet for **MR226675L3V**

### **Rbm41 (NM\_001172147) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Rbm41 (NM_001172147) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rbm41
Synonyms:	BB131279; D330023I21Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001172147
ORF Size:	1275 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226675).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001172147.1</a> , <a href="#">NP_001165618.1</a>
RefSeq Size:	5621 bp
RefSeq ORF:	1278 bp
Locus ID:	237073
UniProt ID:	<a href="#">Q8JZV4</a>
Cytogenetics:	X F1
Gene Summary:	May bind RNA.[UniProtKB/Swiss-Prot Function]



[View online »](#)