

Product datasheet for **MR226681L3V**

Trim8 (NM_053100) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Trim8 (NM_053100) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Trim8
Synonyms:	AA408830; BB006755; GERP; Rnf27
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053100
ORF Size:	1653 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226681).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053100.2 , NP_444330.2
RefSeq Size:	3262 bp
RefSeq ORF:	1656 bp
Locus ID:	93679
UniProt ID:	Q99PJ2
Cytogenetics:	19 C3
Gene Summary:	Probable E3 ubiquitin-protein ligase which may promote proteasomal degradation of SOCS1. [UniProtKB/Swiss-Prot Function]



[View online »](#)