

Product datasheet for **MR226687L3V**

Rbm10 (NM_001167775) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rbm10 (NM_001167775) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rbm10
Synonyms:	E430039K10Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001167775
ORF Size:	2787 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226687).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001167775.1 , NP_001161247.1
RefSeq Size:	3715 bp
RefSeq ORF:	2790 bp
Locus ID:	236732
UniProt ID:	Q99KG3
Cytogenetics:	X A1.3
Gene Summary:	Not known. Binds to RNA homopolymers, with a preference for poly(G) and poly(U) and little for poly(A) (By similarity). May bind to specific miRNA hairpins (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)