

Product datasheet for **MR226694**

Rps6ka3 (NM_148945) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rps6ka3 (NM_148945) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Rps6ka3
Synonyms:	MAPKAPK-1b; MPK-9; p90RSK3; pp90RSK2; Rsk2; S6K-alpha3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR226694 representing NM_148945
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGCCGCTGGCGCAGCTGGCGGACCCGTGGCAGAAGATGGCTGTGGAGAGCCCTCCGACAGCGCGGAGA
ATGGACAGCAAATTATGGATGAACCTATGGGAGAGGAGGAGATTAACCCACAAACTGAAGAAGGCAGTAT
CAAAGAAATTGCAATCACACATCATGTGAAGGAAGGACATGAAAAGGCAGATCCTCCAGTTTGAACCT
TAAAAAGTATTAGGGCAGGGATCATTGGAAAAGTTTTCTTAGTTAAAAAATCTCAGGCTCTGATGCTA
GACAGCTTTATGCCATGAAAGTATTAAGAAGGCCACGCTGAAAGTTCGAGACCGTGTCCGACAAAAAT
GGAACGTGATATCTTGGTAGAAGTCAATCACCTTTTCATTGTCAAATGCATTACGCTTTTCAAACGGAA
GGGAAGTTGTATCTTATTTGGATTTTCTCAGGGCGGAGACTTGTTCACGCTTATCCAAAGAGGTGA
TGTTACAGAGGAAGATGTCAAATCTACTTGGCTGAACCTGCATTGCTTTAGACCATCTTCATAGCCT
GGGAATAATCTATAGAGACTTAAAACCAGAAAACATACTTCTTGATGAAGAAGGTACATCAAGTAACT
GATTTTGGCTTAAGTAAGGAATCTATTGATCATGAGAAGAAGGCTTATTCTTTTTGTGGCACTGTGGAAT
ACATGGCTCCAGAAGTAGTTAACCCGAGAGGTCACACTCAGAGTGCAGGACTGGTGGTCTTTGGTGTGT
GATGTTTGAATGCTCACTGGTACACTACCTTTCCAAGGAAAAGATCGTAAAGAAACAATGACTATGATT
CTTAAAGCCAACTCGGGATGCCACAGTTTCTGAGTCCTGAAGCCAGAGTCTTTTACGAATGCTTTTCA
AACGGAATCCTGCAAACAGATTAGGTGCTGGACCAGATGGAGTTGAAGAAATTAAGACATTCATTTTT
CTCAACAATAGACTGGAATAAATATATAAGAAGAGAGATTCACCCACCTTTTAAAGCCTGCAACTGGCAGA
CCTGAAGATACATTTTATTTGATCCTGAGTTTACTGCAAAAACCTCCAAAGATTCACCGGCATTCCAC
CTAGTGCTAACGCACATCAGCTTTTTCGGGGTTTGTGTTTGTGCTATTACCTCAGATGATGAAAGCCA
AGCTATGCAGACAGTTGGTGTGCATTCAATTGTTGAGCAATTACACAGAACAGTATTGATTTACTGAT
GGATATGAAGTAAAAGAGGATATTGGCGTTGGCTCATACTCCGTTTGTAAAGAGATGTATACATAAGCTA
CAAACATGGAGTTTGGCGTGAAGATTATTGATAAAAGCAAGAGAGACCCAACAGAAGAGATTGAAATCT
TCTTCGCTATGGACAGCATCCAAACATCATTACCCTAAAGGATGTGTATGATGATGGAAAAATGTGTAT
GTAGTAACAGAACTTATGAAAGGAGGTGAATTGCTGGATAAGATTCTTAGACAGAAGTTTTTCTCAGAGC
GAGAGGCCAGCGCTGCTGTTTACTATAACCAAACTGTTGAGTATCTTCATGCACAAGGGTGGTTCA
CAGAGACTTGAACCTAGCAACATCTTTATGTGGATGAGTCTGGTAATCCAGAATCTATTCGAATTTGT
GATTTTGGCTTTGCAAAACAACCTGAGAGCAGAAAATGGTCTTCTCATGACTCCTTGTATACTGCAAAT
TTGTTGCCACAGAGGTTTTAAAACGGCAAGGTTATGATGCTGCCTGTGATATATGGAGTCTTGGCGTCT
CCTTTATACAATGCTTACTGGTTACACTCCATTTGCAAAATGGCCCTGATGATACTCCAGAGGAAACTG
GCACGAATAGGTAGTGGAAAATCTCACTCAGTGGTGGTACTGGAATCTGTTTCAGACACAGCAAAGG
ACCTGGTGTCAAAGATGCTTCATGTAGATCCTCATCAGAGACTGACGGCTGCTCTGGTGTCTCAGACATCC
TTGGATTGTCCACTGGGACCACTACCACAATACCAACTAAACAGACAGGATGCGCCGATCTCGTAAAG
GGTGCCATGGCAGCTACGTACTGCTTTAAACCGCAATCAGTCCCAGTCTTGAACAGTGGGCCGCT
CCACTCTTGCTCAGCGGAGAGGGATTAATAAATCACCTCAACAGCCCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR226694 representing NM_148945
Red=Cloning site Green=Tags(s)


```
MPLAQLADPWQKMAVESPSDSAENGQQIMDEPMGEEEEINPQTEEGSIKEIAITHHVKEGHEKADPSQFEL
LKVLGQGSFGKVFLVKKISGSDARQLYAMKVLKKATLKVRDRVRTKMERDILVEVNHPFIVKLHYAFQTE
GKLYLILDFLRGGDLFTRL SKEVMFTEEDVKFYLAELALALDHLHSLGIIYRDLKPENILLDEEGHIKLT
DFGLSKESIDHEKKAYSFCGTVEYMAPEVVNRRGHTQSADWWSFGVLMFEMLTGTLPFQGKDRKETMTMI
LKAKLGMPQFLSPEAQSLRMLFKRNPANRLGAGPDGVEEIKRHSFFSTIDWNKLYRREIHPPFKPATGR
PEDTFYDFPEFTAKTPKDSPIPPSANAHQLFRGFSVAITSDDESQAMQTVGVHSIVQQLHRNSIQFTD
GYEVKEDI GVGSYVCKRCIHKATNMEFAVKIIDKSKRDPTEEIEILLRYGQHPNIIITLKDVYDDGKYVY
VVELMKGSELLDKILRQKFFSEREASAVLFTITKTVEYLHAQGVVHRDLKPSNIIYVDESGNPESIRIC
DFGFAKQLRAENGLLMTPCYTANFVAPEVLKRQGYDAACDIWSLVLLYMLTGYTPFANGPDDTPEEIL
ARIGSGKFSLSGGYWNSVSDTAKDLVSKMLHVDPHQRLTAALVLRHPWIVHWDQLPQYQLNRQDAPHLVK
GAMAATYSALNRNQSPVLEPVGRSTLAQRRGIKKITSTAL
```

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9009_g03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



* The last codon before the Stop codon of the ORF

ACCN: NM_148945

ORF Size: 2220 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)


OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_148945.2</u> , <u>NP_683747.1</u>
RefSeq Size:	2223 bp
RefSeq ORF:	2223 bp
Locus ID:	110651
UniProt ID:	<u>P18654</u>
Cytogenetics:	X 73.27 cM
MW:	84.1 kDa

Gene Summary:

Serine/threonine-protein kinase that acts downstream of ERK (MAPK1/ERK2 and MAPK3/ERK1) signaling and mediates mitogenic and stress-induced activation of the transcription factors CREB1, ETV1/ER81 and NR4A1/NUR77, regulates translation through RPS6 and EIF4B phosphorylation, and mediates cellular proliferation, survival, and differentiation by modulating mTOR signaling and repressing pro-apoptotic function of BAD and DAPK1. In fibroblast, is required for EGF-stimulated phosphorylation of CREB1 and histone H3 at 'Ser-10', which results in the subsequent transcriptional activation of several immediate-early genes. In response to mitogenic stimulation (EGF and PMA), phosphorylates and activates NR4A1/NUR77 and ETV1/ER81 transcription factors and the cofactor CREBBP. Upon insulin-derived signal, acts indirectly on the transcription regulation of several genes by phosphorylating GSK3B at 'Ser-9' and inhibiting its activity. Phosphorylates RPS6 in response to serum or EGF via an mTOR-independent mechanism and promotes translation initiation by facilitating assembly of the preinitiation complex. In response to insulin, phosphorylates EIF4B, enhancing EIF4B affinity for the EIF3 complex and stimulating cap-dependent translation. Is involved in the mTOR nutrient-sensing pathway by directly phosphorylating TSC2 at 'Ser-1798', which potently inhibits TSC2 ability to suppress mTOR signaling, and mediates phosphorylation of RPTOR, which regulates mTORC1 activity and may promote rapamycin-sensitive signaling independently of the PI3K/AKT pathway. Mediates cell survival by phosphorylating the pro-apoptotic proteins BAD and DAPK1 and suppressing their pro-apoptotic function. Promotes the survival of hepatic stellate cells by phosphorylating CEBPB in response to the hepatotoxin carbon tetrachloride (CCl4). Is involved in cell cycle regulation by phosphorylating the CDK inhibitor CDKN1B, which promotes CDKN1B association with 14-3-3 proteins and prevents its translocation to the nucleus and inhibition of G1 progression. In LPS-stimulated dendritic cells, is involved in TLR4-induced macropinocytosis, and in myeloma cells, acts as effector of FGFR3-mediated transformation signaling, after direct phosphorylation at Tyr-529 by FGFR3. Phosphorylates DAPK1 (By similarity). Negatively regulates EGF-induced MAPK1/3 phosphorylation via phosphorylation of SOS1. Phosphorylates SOS1 at 'Ser-1134' and 'Ser-1161' that create YWHAB and YWHAЕ binding sites and which contribute to the negative regulation of MAPK1/3 phosphorylation (PubMed:22827337). Phosphorylates EPHA2 at 'Ser-897', the RPS6KA-EPHA2 signaling pathway controls cell migration (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR226694