

## Product datasheet for **MR226695L1V**

### Sh2d1a (NM\_011364) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Sh2d1a (NM_011364) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sh2d1a
Synonyms:	Gm686; SAP
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_011364
ORF Size:	378 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226695).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_011364.3</a> , <a href="#">NP_035494.1</a>
RefSeq Size:	819 bp
RefSeq ORF:	381 bp
Locus ID:	20400
UniProt ID:	<a href="#">O88890</a>
Cytogenetics:	X A4



[View online »](#)

**Gene Summary:**

Cytoplasmic adapter regulating receptors of the signaling lymphocytic activation molecule (SLAM) family such as SLAMF1, CD244, LY9, CD84, SLAMF6 and SLAMF7. In SLAM signaling seems to cooperate with SH2D1B/EAT-2. Initially it has been proposed that association with SLAMF1 prevents SLAMF1 binding to inhibitory effectors including INPP5D/SHIP1 and PTPN11/SHP-2. However, by simultaneous interactions, recruits FYN which subsequently phosphorylates and activates SLAMF1 (By similarity). Positively regulates CD244/2B4- and CD84-mediated natural killer (NK) cell functions (PubMed:22683124). Can also promote CD48-, SLAMF6-, LY9-, and SLAMF7-mediated NK cell activation (PubMed:19648922). In the context of NK cell-mediated cytotoxicity enhances conjugate formation with target cells (PubMed:22683124). May also regulate the activity of the neurotrophin receptors NTRK1, NTRK2 and NTRK3.[UniProtKB/Swiss-Prot Function]