

Product datasheet for **MR226717L3V**

Nfe2l2 (NM_010902) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nfe2l2 (NM_010902) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nfe2l2
Synonyms:	Nr; Nrf2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010902
ORF Size:	1791 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226717).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010902.3
RefSeq Size:	2469 bp
RefSeq ORF:	1794 bp
Locus ID:	18024
UniProt ID:	Q60795
Cytogenetics:	2 44.75 cM



[View online »](#)

Gene Summary:

This gene encodes a transcription factor which is a member of a small family of basic leucine zipper (bZIP) proteins. The encoded transcription factor regulates genes which contain antioxidant response elements (ARE) in their promoters; many of these genes encode proteins involved in response to injury and inflammation which includes the production of free radicals. [provided by RefSeq, Sep 2015]