

Product datasheet for **MR226733L4V**

Gfap (NM_001131020) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gfap (NM_001131020) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gfap
Synonyms:	AI836096
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001131020
ORF Size:	1284 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226733).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001131020.1 , NP_001124492.1
RefSeq Size:	2600 bp
RefSeq ORF:	1287 bp
Locus ID:	14580
UniProt ID:	P03995
Cytogenetics:	11 66.48 cM
Gene Summary:	GFAP, a class-III intermediate filament, is a cell-specific marker that, during the development of the central nervous system, distinguishes astrocytes from other glial cells. [UniProtKB/Swiss-Prot Function]



[View online »](#)