

## Product datasheet for **MR226742L3V**

### **Kcna3 (NM\_008418) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Kcna3 (NM_008418) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Kcna3
Synonyms:	Kca1-3; Kv1.3; Mk-3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008418
ORF Size:	1587 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226742).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_008418.2</a> , <a href="#">NP_032444.2</a>
RefSeq Size:	1968 bp
RefSeq ORF:	1587 bp
Locus ID:	16491
UniProt ID:	<a href="#">P16390</a>
Cytogenetics:	3 46.59 cM



[View online »](#)

**Gene Summary:**

Mediates the voltage-dependent potassium ion permeability of excitable membranes. Assuming opened or closed conformations in response to the voltage difference across the membrane, the protein forms a potassium-selective channel through which potassium ions may pass in accordance with their electrochemical gradient.[UniProtKB/Swiss-Prot Function]