

Product datasheet for **MR226743A1V**

Mouse Etv1 (NM_001163154) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Etv1 (NM_001163154) AAV Particle
Tag:	Myc-DDK
Symbol:	Etv1
Synonyms:	ER81; Etsrp81
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR226743 representing NM_001163154 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCTTCAAGATTTAAGTGCAGGCGTCTTCTCCCTCCTTGTTCAACACAGAACATTAGCTCAGGTAC
CTGACAATGATGAGCAGTTTGTTCAGACTATCAGGCTGAAAGTTTGGCTTTTCATGGCCTCCCACTGAA
AATCAAGAAAGAACCCACAGCCCATGTTCAAGAACTCGGGTCTGCTTGCAGTCAAGAGCAGCCCTTTAA
TTCAGCTATGGAGAAAAGTGCTGTACAATGTCAGTGCCTATGATCAGAAGCCACAAGTGGGAATGAGGC
CCTCCAATCCCCGACACCATCCAGCACTCCAGTGTCCCACATACACCATGCATCTCCAAACACAGCTCA
TACTCCGAAACCTGACCGGGCCTTCCCAGCTCACCTCCCTCCTTCTCAGTCCATACCAGACAGCACCTAC
CCCATGGACCACAGATTTGCGCCGCCAGCTTTCTGAACCCTGTAATTCTTTTCTCCTTTGCCGACAATGC
CAAGGGAAGGGCGTCTATGTACCAACGCCAGATGTCTGAGCCAAACATCCCCTTCCCACCACAAGGCTT
TAAGCAGGAGTACCACGACCCAGTCTATGAACATACCACCATGGTTGGCGGTGCAGCCAGCCAAAGCTTC
CCCCCTCCTTTGATGATCAAACAGGAACCCAGAGATTTTGCATATGACTCAGAAGTGCCTAGCTGCCACT
CCATTTATATGAGGCAAGAAGGCTTCTGGCTCATCAAGCAGAACAGAAGGCTGCATGTTTGAGAAGGG
TCCCAGGCAGTTCTATGATGACACCTGTGTGGTCCCGAGAAGTTCGATGGGGACATTAAGCAAGAGCCT
GGAATGTACCGGAAGGACCCACGTACCAGAGGCGAGGATCCCTCCAGCTCTGGCAGTTTTTGGTAGCTC
TTCTGGATGACCTTCAAATTCTCATTTTCATTGCTGGACTGGACGAGGCATGGAATTTAACTGATTGA
GCCCCAAGAGGTGGCCCGGCTTGGGGCATTGAGAAGAAGGCGGCGATGAACATGACAAACTTAGT
CGTTCTCTCCGCTATTATTATGAGAAGGGAATCATGAAAAGGTGGCTGGAGAAAAGATACGTGTACAAAT
TTGTGTGTGACCCGGAAGCCCTTTTCTATGGCCTTTCGGATAACACGCGCCCACTGCTGAAGACGGA
CATGGAACGTACATCAACGAAGAGGACACAGTGCCTCTGTCTCACTTTGATGAGAGCATGACCTACATG
CCCGAAGGGGCTGCTGCAACCTCACCCCTACAACGAAGGATACGTGTAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR226743 representing NM_001163154 Red=Cloning site Green=Tags(s) MLQDLSAGVFFPPCSQHRTLAQVPDNEQFVPDYQAESLAFHGLPLKIKKEPHSPCSELGSACSQEPPFK FSYGEKCLYNVSAYDQKPQVGMRPSNPPTPSSTPVSPHASPNTAHTPKPDRAFPAPHLPPSQSIPDSTY PMDHRFRRLSEPCNSFPPLPTMPREGRPYQRMSEPNIFFPPQGFQYHDPVYEHTTMVGGAASQSF PPPLMIKQEPDFAYDSEVPSCHSIYMRQEGFLAHPSTRTEGCMFEKGPRQFYDDTCVVPEKFDGDIKQEP GMYREGPTYQRRGSLQLWQFLVALLDDPSNSHFIAWTGRGMEFKLIEPEEVARRWGIQKNRPAMNYDKLS RSLRYYYEKGIMQKVAGERYVYKFVCDPEALFSMAFPDNQRPLLKTDMERHINEEDTVPLSHFDESMTYM PEGGCCNPHYNEGYVY TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_001163154
ORF Size:	1311 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001163154.1</u> , <u>NP_001156626.1</u>
RefSeq Size:	3806 bp
RefSeq ORF:	1314 bp
Locus ID:	14009
Cytogenetics:	12 17.99 cM
MW:	50.6 kDa