

## Product datasheet for MR226748L4

### Il4 (NM\_021283) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Il4 (NM_021283) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Il4
Synonyms:	BSF-1; Il-4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226748).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.


ACCN:	NM_021283
ORF Size:	420 bp



[View online »](#)

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_021283.2</a> , <a href="#">NP_067258.1</a>
<b>RefSeq Size:</b>	605 bp
<b>RefSeq ORF:</b>	423 bp
<b>Locus ID:</b>	16189
<b>UniProt ID:</b>	<a href="#">P07750</a>
<b>Cytogenetics:</b>	11 31.97 cM
<b>Gene Summary:</b>	Participates in at least several B-cell activation processes as well as of other cell types (PubMed:3083412). It is a costimulator of DNA-synthesis. It induces the expression of class II MHC molecules on resting B-cells (PubMed:3498301). It enhances both secretion and cell surface expression of IgE and IgG1 (PubMed:3498301). It also regulates the expression of the low affinity Fc receptor for IgE (CD23) on both lymphocytes and monocytes. Positively regulates IL31RA expression in macrophages (PubMed:25847241). Stimulates autophagy in dendritic cells by interfering with mTORC1 signaling and through the induction of RUFY4 (PubMed:26416964).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR226748L4