

Product datasheet for **MR226766L4V**

P2ry1 (NM_008772) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	P2ry1 (NM_008772) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	P2ry1
Synonyms:	P2Y1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008772
ORF Size:	1122 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226766).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008772.4 , NP_032798.1
RefSeq Size:	3899 bp
RefSeq ORF:	1122 bp
Locus ID:	18441
UniProt ID:	P49650
Cytogenetics:	3 E1
Gene Summary:	Receptor for extracellular adenine nucleotides such as ADP. In platelets, binding to ADP leads to mobilization of intracellular calcium ions via activation of phospholipase C, a change in platelet shape, and ultimately platelet aggregation.[UniProtKB/Swiss-Prot Function]



[View online »](#)