

Product datasheet for **MR226827L3V**

Csf1 (NM_007778) Mouse Tagged ORF Clone Lentiviral Particle

Product data:


Product Type:	Lentiviral Particles
Product Name:	Csf1 (NM_007778) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Csf1
Synonyms:	C87615; Csfm; MCSF; op
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007778
ORF Size:	1656 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226827).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007778.4 , NP_031804.3
RefSeq Size:	4192 bp
RefSeq ORF:	1659 bp
Locus ID:	12977
UniProt ID:	P07141
Cytogenetics:	3 46.83 cM



[View online »](#)

Gene Summary:

Cytokine that plays an essential role in the regulation of survival, proliferation and differentiation of hematopoietic precursor cells, especially mononuclear phagocytes, such as macrophages and monocytes. Promotes the release of proinflammatory chemokines, and thereby plays an important role in innate immunity and in inflammatory processes. Plays an important role in the regulation of osteoclast proliferation and differentiation, the regulation of bone resorption, and is required for normal bone development. Required for normal male and female fertility. Promotes reorganization of the actin cytoskeleton, regulates formation of membrane ruffles, cell adhesion and cell migration. Plays a role in lipoprotein clearance. [UniProtKB/Swiss-Prot Function]

Product images:

[MR226827L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with MR226827L3V particle to overexpress human Csf1-Myc-DDK fusion protein.