

Product datasheet for **MR226828L3V**

Nfe2l1 (NM_008686) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nfe2l1 (NM_008686) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nfe2l1
Synonyms:	AA408798; AW212678; LCR-F1; Lcrf1; NRF1; TCF-11; TCF11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008686
ORF Size:	2223 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226828).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008686.3 , NP_032712.2
RefSeq Size:	4654 bp
RefSeq ORF:	2226 bp
Locus ID:	18023
UniProt ID:	Q61985
Cytogenetics:	11 60.12 cM



[View online »](#)

Gene Summary:

This gene encodes a protein that is involved in globin gene expression in erythrocytes. Confusion has occurred in bibliographic databases from the use of NRF1 for this gene, NFE2L1, and for "nuclear respiratory factor 1" which has an official symbol of NRF1. [provided by RefSeq, Jul 2008]