

## Product datasheet for **MR226847L3V**

### Vav1 (NM\_001163815) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Vav1 (NM_001163815) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Vav1
Synonyms:	Vav; vav-T
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001163815
ORF Size:	2463 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226847).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001163815.1</a> , <a href="#">NP_001157287.1</a>
RefSeq Size:	4096 bp
RefSeq ORF:	2466 bp
Locus ID:	22324
Cytogenetics:	17 29.76 cM
Gene Summary:	Couples tyrosine kinase signals with the activation of the Rho/Rac GTPases, thus leading to cell differentiation and/or proliferation.[UniProtKB/Swiss-Prot Function]



[View online »](#)