

Product datasheet for **MR226848L3V**

Bdkrb1 (NM_007539) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Bdkrb1 (NM_007539) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Bdkrb1
Synonyms:	B1BKR; B1R; Bdkrb; BKR1; BRADYB1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007539
ORF Size:	1005 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226848).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007539.2
RefSeq Size:	1338 bp
RefSeq ORF:	1005 bp
Locus ID:	12061
UniProt ID:	Q61125
Cytogenetics:	12 E
Gene Summary:	This is a receptor for bradykinin. Could be a factor in chronic pain and inflammation. [UniProtKB/Swiss-Prot Function]



[View online »](#)